

- 14GA 2 x 2" STEEL TUBE CONSTRUCTION TYPICAL
- ASSEMBLY REQUIRED
- STATIC LOAD CAPACITY, EQUALLY DISTRIBUTED
  - COMBINED SYSTEM, NO CASTERS : 1200 LBS (880 LBS W/ CASTERS)
  - LOWER SHELF : 200 LBS
  - UPPER SHELVES : 150 LBS
- AVAILABLE SEPARATELY
  - UPRIGHT POST KIT WITH HORIZONTAL SUPPORTS (SHOWN)
  - 12" DEEP SHELF KIT (SHOWN WITH 2 SHELF KITS)
  - CHEM RESISTANT PHENOLIC WORK SURFACE (SHOWN)
  - SIDE APRON TO ACCOMMODATE 36" DEEP CUSTOM WORK SURFACE
  - SINGLE SUSPENDED DRAWER
  - POWER STRIP KIT (REQ'S UPRIGHT POST KIT)
  - BIN BOX MOUNT BRACKET (REQ'S UPRIGHT POST KIT)
  - SWIVEL CASTER KIT W/BRAKE
  - 24.75" DEEP LOWER SHELF W/BRACKETS (SHOWN)
  - SAFETY DIVIDER KITS; SIDE AND REAR, 1/4" ACRYLIC
- SEE PRICE LIST FOR FRAME AND ACCESSORY PART NUMBERS



K13409-48: 48W UPRIGHT POST KIT  
 K13409-72: 72W UPRIGHT POST KIT  
 K13032-12: 48x12 UPPER SHELF W/ BKTS  
 K13034-12: 72x12 UPPER SHELF W/ BKTS



LL4830CRP: 48x30 PHENOLIC TOP  
 LL7230CRP: 72x30 PHENOLIC TOP

K7003-48: 48W LOWER SHELF W/ BKTS  
 K7003-72: 72W LOWER SHELF W/ BKTS

UTFH48B: FRAME, 48W FIXED HT  
 UTFH72B: FRAME, 72W FIXED HT

STANDARD SIZE MATRIX		
Table size	Frame width	Worksurface width
48 x 30	47-3/4"	48"
72 x 30	71-3/4"	72"

DRAWN	4/20/2020	LOOPED LOGIC		
CHECKED		747 Calle Plano, Camarillo CA 93012   866-340-2622		
QA		TITLE CUT SHEET, LOOPED LOGIC UTILITY TABLE, FIXED HT		
APPROVED		FINISH PETAL WHITE		
UNLESS OTHERWISE NOTED, GEN. TOLERANCE +/- .015"; ALL DIMENSIONS ARE IN INCHES	JOB #	SIZE	PART NO	REV
	QUOTE #	D		A
		SCALE	.15 MATL	SHEET 1 OF 1

