






## *AGILITY™ SERVICE HUB modules*




Type 1




Type 2



Type 3



Type 4



Type 5

*Unique power, data and specialty gas service utility modules*

*Designed and manufactured in Camarillo, California* 

[www.HANSONLAB.com](http://www.HANSONLAB.com) | Camarillo, CA | 805-498-3121

## Features:

- Choices for AC power, RJ45 ethernet, specialty gases, USB power, isolated/multiple circuits and dedicated emergency power
- Can be tailored for your specific needs and pre-installed into your new-build HLS table frame, resulting in a clean, uncompromised work center
- Exclusive design distributes power and services vertically, for the most efficient use of the workspace
- Extruded aluminum-alloy housing with clear anodized finish delivers superior durability and aesthetics
- Certified to UL and CSA standards



## Design for your application:

- To configure an HLS Service Hub (SH) module for your specific requirements, begin by selecting a faceplate type (visit the [HLS website](#) to see the available options, and download a cut-sheet).
- Next, configure the component details of your SH module by referencing the table on the HLS cut-sheet for your specific faceplate type and configuration.
- Finally, identify your specific SH part number based on the part numbering key at the bottom of the cut-sheet reference page.

