



AGILITY™ ADAPTABLE workstation



Heavy-duty adjustable height bench systems

Designed and manufactured in Camarillo, California 

Features:

- Reconfigurable modular design features adjustable worktop height
- Heavy-gauge welded steel frame supports 3500+ lb. working capacity
- Integrated variable-height shelving system for flexibility
- Choose from multiple lengths, heights and depths
- Available Service Hub system -- offers numerous power, data and specialty gas services - can be tailored for your specific needs and pre-installed into the frame, resulting in a clean, uncompromised work center

Build it:

To configure an HLS Agility™ Adaptable bench for your specific requirements, begin by making a few selections:

- Workbench length
- Front-to-rear depth
- Upright height
- Choose your optional Service Hub(s)
 - Unique power, data and specialty gas service utility modules
 - 18 different configurations using five distinct types; provides numerous service combinations
 - Refer to separate HLS Service Hub (SH) literature for details on the range of available options

Agility™ Adaptable

Table system part numbering guide

AA- 60 30 A 84 R

AGILITY ADAPTABLE]	AA- 60 30 A 84 R
Work top length]	60 - 60"
	36 - 36"
	48 - 48"
	72 - 72"
Work top depth]	30 - 30"
	36 - 36"
	42 - 42"
	60 - 60" (double sided table)
Work top height]	A - Adjustable (3 position; 30.75" - 36.75")
Upright height]	84 - 84"
	78 - 78"
	72 - 72"
Service Hub location]	R - Right post mounting
	L - Left post mounting
	B - Both; one SH mounting each post
	N - None. SH delete

P/N example for 60"W x 30"D x 84"H bench with right-hand Service Hub: **AA-6030A84R**

