



AGILITY™ ESSENTIAL workstation



Heavy-duty modular bench systems

Designed and manufactured in Camarillo, California 🇺🇸

Features:

- Heavy-gauge welded steel frame supports 3500+ lb. working capacity
- Configurable, modular design supports future expansion
- Integrated adjustable-height shelving system for flexibility
- Multiple lengths, heights and depths available to meet your specifications
- Available Service Hub system -- offers numerous power, data and specialty gas services - can be tailored for your specific needs and pre-installed into the frame, resulting in a clean, uncompromised work center

Build it:


To configure an HLS Agility™ Essential bench for your specific requirements, begin by making a few selections:

- Workbench length
- Front-to-rear depth
- Worktop height
- Upright height
- Choose your optional Service Hub(s)
 - Unique power, data and specialty gas service utility modules
 - 18 different configurations using five distinct types; provides numerous service combinations
 - Refer to separate HLS Service Hub (SH) literature for details on the range of available options

Agility™ Essential

Table system part numbering guide

AE- 48 30 H 84 B



P/N example for 48"W x 30"D x 84"H bench
 with two Service Hubs: **AE-4830H84B**

