



## *AGILITY™ MOBILE workstation*



*Portable heavy-duty bench systems*

*Designed and manufactured in Camarillo, California* 

### Features:

- Thick-gauge welded steel frame supports 3500+ lb. working capacity
- Heavy-duty locking casters and adjustable leveling feet provide superior stability
- Integrated adjustable-height shelving system for flexibility
- Configurable in multiple lengths, heights and depths
- Available Service Hub system -- offers numerous power, data and specialty gas services - can be tailored for your specific needs and pre-installed into the frame, resulting in a clean, uncompromised work center

**Build it:**

To configure an HLS Agility™ Mobile bench for your specific requirements, begin by making a few selections:

- Workbench length
- Front-to-rear depth
- Upright height
- Fixed or adjustable worktop height
- Choose your optional Service Hub(s)
  - Unique power, data and specialty gas service utility modules
  - 18 different configurations using five distinct types; provides numerous service combinations
  - Refer to separate HLS Service Hub (SH) literature for details on the range of available options

**Agility™ Mobile**

Table system part numbering guide

**AM- 60 30 H 84 R**

**AGILITY MOBILE ]**

**Work top length ]**

- 36 - 36"
- 48 - 48"
- 60** - 60"
- 72 - 72"
- 96 - 96" (2-bay; 3-upright)

**Work top depth ]**

- 30** - 30"
- 36 - 36"
- 42 - 42"

**Work top height ]**

- H** - High (standing ht); fixed (36.75" nom)
- A - Adjustable (30.75" - 36.75")

**Upright height ]**


- 72 - 72"
- 78 - 78"
- 84** - 84"

N - None. Upright delete (N/A on adjustable height tables)


**Service Hub location ]**

- L - Left post mounting
- R** - Right post mounting
- B - Both; one SH mounting each post
- C - Center post mounting (96" long tables only)
- A - Left, right and center post mounting (96" long tables only)
- N - None. SH delete

P/N example for 60"W x 30"D x 84"H, fixed-height bench with right-hand Service Hub: **AM-6030H84R**



12" DEEP, ADJUSTABLE HEIGHT SHELVES:  
 -72" HIGH MODELS INCLUDE ONE ROW  
 -78", 84" HIGH MODELS INCLUDE TWO ROWS (SHOWN)



OPTIONAL SERVICE HUB (SH) AND POWER CORD SHOWN - SEE SH PLANNING GUIDE FOR DETAILS



HEIGHT ADJUSTABLE WORKTOP CONFIGURATION AVAILABLE

NOTE: 96" LONG TABLES ARE CONFIGURED WITH 3 UPRIGHTS AND 2 ROWS OF 4 SHELVES