

- 14GA 2 x 2" STEEL TUBE CONSTRUCTION TYPICAL
- ASSEMBLY REQUIRED
- SEE PRICE LIST FOR ADDITIONAL PART NUMBERS
- AVAILABLE SEPARATELY-
  - UPRIGHT POST KIT WITH HORIZONTAL SUPPORTS (SHOWN)
  - 12" DEEP SHELF KIT (SHOWN WITH 2 SHELF KITS)
  - CHEM RESISTANT PHENOLIC WORKSURFACE, 3/4" (SHOWN)
  - SIDE APRON TO ACCOMMODATE 36" DEEP CUSTOM WORKSURFACE
  - SINGLE SUSPENDED DRAWER
  - 3" DIA STEM SWIVEL CASTER KIT W/BRAKE
  - 24.75" DEEP LOWER SHELF W/BRACKETS (SHOWN)
  - SAFETY DIVIDER KITS; SIDE AND REAR, 1/4" ACRYLIC



- K13409-48: 48W UPRIGHT POST KIT
- K13409-72: 72W UPRIGHT POST KIT
- K13032-12: 48x12 UPPER SHELF W/ BKTS
- K13034-12: 72x12 UPPER SHELF W/ BKTS



- LL4830CRP: 48x30 PHENOLIC TOP
- LL7230CRP: 72x30 PHENOLIC TOP

- K7003-48: 48W LOWER SHELF W/ BKTS
- K7003-72: 72W LOWER SHELF W/ BKTS

- UTFH48B: FRAME, 48W FIXED HT
- UTFH72B: FRAME, 72W FIXED HT

STANDARD SIZE MATRIX		
Table size	Frame width	Worksurface width
48 x 30	47-3/4"	48"
72 x 30	71-3/4"	72"

DRAWN	MWR	4/20/2020	LOOPED LOGIC		
CHECKED			747 Calle Plano, Camarillo CA 93012   866-340-2622		
QA			TITLE		
			CUT SHEET, LOOPED LOGIC UTILITY TABLE, FIXED HT		
			FINISH		
			PETAL WHITE		
APPROVED			SIZE	PART NO	REV
			D		A
UNLESS OTHERWISE NOTED, GEN. TOLERANCE +/- .015"; ALL DIMENSIONS ARE IN INCHES		JOB #	SCALE	.15	MATL
		QUOTE #			SHEET 1 OF 1

