

AGILITY™ ESSENTIAL WORKSTATION



Introductory Line

Improve your bottom line with expandable workstations.

FEATURES:

- Heavy-gauge welded steel frame supports 3500+ lb. working capacity
- Configurable, modular design supports future expansion
- Integrated adjustable-height shelving system for flexibility
- Multiple lengths, heights and depths available to meet your specifications
- Tables can be linked together for a professional touch and expanded working space

BUILD YOURS:

Select AE when choosing your desk specifications. Make sure to also choose your length, depth, height, and Service Hub Locations. Available options to the right.

AE- _____

Agility™ Essential - stand alone
Table system part numbering guide

AE- 48 30 H 84 B

AGILITY ESSENTIAL]

Work top length]
 48 - 48"
 60 - 60"
 72 - 72"

Work top depth]
 30 - 30"
 36 - 36"
 60 - 60" (double sided table)

Work top height]
 L - Low (sitting ht); fixed (30.75" nom)
 M - Medium (ADA ht); fixed (32.75" nom)
 H - High (standing ht); fixed (36.75" nom)

Upright height]
 78 - 78"
 84 - 84"
 N - None. Upright delete

Service Hub location]
 L - Left post mounting
 R - Right post mounting
 B - Both; one SH mounting each post
 N - None. SH delete

P/N example for 48"W x 30"D x 84"H bench
with two Service Hubs: **AE-4830H84B**

