



Hanson M3 & M4 Series Workstations

The newest members of the HLS M-Series

HLS M3 and M4 bench systems are highly configurable, and offer numerous accessory options designed to meet a variety of project budgets.

The backbone of the M3 & M4 bench is a robust 2"x 5" 11-gauge powder-coated steel upright for maximum stability and extended life.


Sizing choices:

- Widths – 4', 5' and 6'
- Depths - 30", 36", 60" (double-sided)
- Height adjustable in 1" increments from 30-3/4 to 36-3/4"

Power, data, and gas service options can be pre-installed into the frame, resulting in a clean, uncompromised work center.

Additional accessories:



- 1" thick Lab-grade phenolic work surfaces
- Shelving - multiple depth choices with various retaining lip styles
- Mobile and suspended cabinets - numerous door and drawer configurations
- Bin box mount
- Horizontal power strip



M3 bench (single-sided)



M4 bench (double-sided)



Designed and manufactured in Camarillo, California

Hanson Lab Solutions - 805.498.3121

www.HANSONLAB.com

