

# AGILITY™ ADAPTABLE WORKSTATION



**Fully Adjustable Countertop** - Elevate your work with a table that is just as versatile as you are.

## **FEATURES:**

- Reconfigurable modular design features adjustable worktop height
- Heavy-gauge welded steel frame supports 3500+ lb. working capacity
- Integrated variable-height shelving system for flexibility
- 6-position counter top (30" to 36" Adjustable)

**BUILD YOURS:**

Select AA when choosing your desk specifications. Make sure to also choose your length, depth, height, and Service Hub Locations. Available options to the right.

**AA-** \_\_\_\_\_

*Agility™ Adaptable*  
Table system part numbering guide

**AA- 60 30 A 84 R**

**AGILITY ADAPTABLE ]**

**Work top length ]**  
36 - 36"  
48 - 48"  
**60 - 60"**  
72 - 72"

**Work top depth ]**  
**30 - 30"**  
36 - 36"  
42 - 42"  
60 - 60" (double sided table)

**Work top height ]**  
**A - Adjustable (6 position; 36.75" - 30.75")**

**Upright height ]**  
72 - 72"  
78 - 78"  
**84 - 84"**

**Service Hub location ]**  
L - Left post mounting  
**R - Right post mounting**  
B - Both; one SH mounting each post  
N - None. SH delete

P/N example for 60"W x 30"D x 84"H bench with right-hand Service Hub: **AA-6030A84R**

