




Ceiling Service Panels

Ceiling Service Panels


- Designed to fit inside a standard ceiling tile grid when service connection from the ceiling is desired
- Typically paired with modular benches equipped with Service Hub or accessory panel
- Painted steel panel construction
- Stainless Steel available (-SS)
- Service panels may be customized to meet your project requirements
- Order service fixtures and electrical / data components separately




C1600
23-1/2" x 23-1/2" 8 Electrical cutouts & 8 Plumbing cutouts




C1601
23-1/2" x 23-1/2" 8 Electrical cutouts & 4 Plumbing cutouts



C1603
23-1/2" x 11-1/2" 4 Electrical cutouts & 4 Plumbing cutouts



C1605
23-1/2" x 47-1/2" 10 Electrical cutouts & 6 Plumbing cutouts



C1606
11-1/2" x 47-1/2" 8 Electrical cutouts & 8 Plumbing cutouts