

Workstations and table options for every need, budget, and design solution.

AGILITY: Our premier bench with max load carrying capacity - features the HLS Service Hub option.

M-SERIES: 2 x 3" or 2 x 5" frame uprights; available with essential power and specialty gas services.

HDW: 2 x 2" welded construction; available with integrated shelving upright posts.

UTILITY: 2 x 2" knockdown design; optional bolt-on shelving upright posts.

HLS shelving and accessories are interchangeable within most HLS table families.



## Adaptable Furniture & Table Index

Page

**WORKSTATION / TABLE COMPARISON CHART** ..... 70

**AGILITY WORKSTATIONS** .....

- Agility Mobile ..... 71
- Agility Adaptable ..... 72
- Agility Essential ..... 73
- Agility Service Hub ..... 75
- Agility Power, Gas and Data Options ..... 78

**M-SERIES WORKSTATIONS** .....

- Features and Options ..... 79
- Numbering Guide..... 80
- Choice Sheets..... 81
- Pre-Configured Bench Assemblies ..... 89
- Accessories..... 90
- M1/M2 Power Packs..... 91
- M3/M4 Power Packs..... 92
- M-Series Power, Gas, and Data Options..... 94

**HEAVY DUTY WELDED TABLES** .....

- Fixed Height ..... 96
- Adjustable Height ..... 101

**UTILITY TABLES** .....

- Utility Table Frames ..... 105
- Utility Table Frame Accessories..... 106

**ACCESSORIES**..... 107

- Table Accessories..... 108
- Ceiling Service Panels and Couplers ..... 110



# ADAPTABLE FURNITURE & TABLE COMPARISON CHART

	Utility	HDW	M-Series	Agility Essential	Agility Adaptable	Agility Mobile
<b>Description</b>	2x2, 14-ga knockdown construction w/ available add-on shelving upright posts	2x2, 14-ga fully-welded construction Available shelving upright posts. Suspended cabinet option.	2x3 (14-ga.) or 2x5 (11-ga.) integrated shelving upright posts. Suspended cabinet option.	2x5, 11-gauge construction w/ integrated shelving upright post available	2x5, 11-gauge construction w/ integrated shelving upright posts. Suspended cabinet option.	2x5, 11-gauge construction w/ integrated shelving upright post available
<b>Static Load Capacity*</b>	1200 lbs.	1200 lbs.	1200 lbs.	3000 lbs.	3000 lbs.	3000 lbs.
<b>Power, Data, Gas Services</b>	Bolt-on horizontal power strip available	Bolt-on horizontal power strip available	Post-integrated pre-wire, pre-plumb available	Full-range pre-wire, pre-plumbed services available w/ HLS Service Hub	Full-range pre-wire, pre-plumbed services available w/ HLS Service Hub	Full-range pre-wire, pre-plumbed services available w/ HLS Service Hub
<b>Mobility</b>	Available caster kit for light-duty applications	Available caster kit for light-duty applications	Modular design allows for expansion	Modular design allows for expansion	Moveable	Fully mobile – move your workstation with instrumentation intact
<b>Countertop Height</b>	Sitting, standing or adjustable height	Sitting, standing or adjustable height	6-position adjustable height	Multiple fixed heights available	6-position adjustable height	Fixed or adjustable height
<b>Widths (nominal)</b>	36", 48", 60", 72"	36", 48", 60", 72", 96"	36", 48", 60", 72"	36", 48", 60", 72"	36", 48", 60", 72", 96"	36", 48", 60", 72", 94"
<b>Depths (nominal)</b>	30, 36"	25", 30", 36", 42"	30", 36", 60" (double-sided)	30", 36", 60" (double-sided)	30", 36", 42", 60" (double-sided)	30", 36", 42"

\*weight evenly distributed, on leveling glides

**Work Surfaces:** Lab-grade phenolic, epoxy resin, plastic laminate, stainless steel (see page 127)

**Shelving:** Multiple depth choices with various retaining lip styles (see page 103)


## Agility Mobile Workstations

Agility Mobile workstations provide maximum flexibility to your room layout. Available with fixed or adjustable height countertop support, you can customize your workstation with a variety of accessories found on the next few pages to create the exact instrument cart or workstation you require.

- Thick-gauge knockdown design (-KD) welded steel frame supports 3000+ lb. capacity
- Heavy-duty locking casters and adjustable leveling feet provide superior stability
- Integrated adjustable-height shelving system for flexibility
- Configurable in multiple lengths, heights and depths
- Countertops heights are adjustable in 1" increments (Available on the adjustable height workstation)
- Assembly required

### Available options:

- Pre-wired / pre-plumbed Service Hub system offers numerous power, data and specialty gas services
- 1" thick lab-grade phenolic work surface (Epoxy or stainless also available)
- Upper shelving - multiples depth choices with various retaining lip styles
- Mobile cabinets - numerous door and drawer configurations
- Horizontal power strip
- Bin box mount



### Agility™ Mobile

Table system part numbering guide

## AM-60 30 H 84 R-KD

- AGILITY MOBILE ]** ———
- Work top length ]** ———
  - 36 - 36"
  - 48 - 48"
  - 60 - 60"
  - 72 - 72"
  - 94 - 93.75" (2-bay; 3-post)
- Work top depth ]** ———
  - 30 - 30"
  - 36 - 36"
  - 42 - 42"
- Work top height ]** ———
  - H - High (standing ht); fixed (36.75" nom)
  - A - Adjustable (36.75" - 30.75")
- Upright post height ]** ———
  - 72 - 72"
  - 78 - 78"
  - 84 - 84"
  - N - None. Post delete (N/A on adj. ht tables)
- Service Hub location\* ]** ———
  - L - Left post mount only
  - R - Right post mount only
  - B - Both; one SH mount each L/R post
  - N - None. SH delete

P/N example for 60"W x 30"D x 84"H, fixed-height bench with right-hand Service Hub:  
AM-6030H84R-KD

\* 93.75" tables only

Use two letter indicators on 93.75" tables to identify Hub locations(s) - one letter per bay. First letter identifies LH bay configuration and second letter identifies RH bay configuration.

P/N example:

93.75"W x 30"D x 84"H fixed-height bench  
LH bay: Service Hub on BOTH post faces  
RH bay: Service Hub on LEFT post face only:  
AM-9430H84BL-KD




# AGILITY ADAPTABLE WORKSTATIONS

## Agility Adaptable Workstations


- Adjustable work top height ( 6 position: 36.75" - 30.75")
- 11ga steel upright posts
- Knockdown design (assembly required)
- Static load capacity, equally distributed combined system: 3000 lbs
- Single or Double-sided configurations

### Available options:

- Pre-wired / pre-plumbed Service Hub system offers numerous power, data and specialty gas services
- 1" thick lab-grade phenolic work surface (Epoxy or stainless also available)
- Upper shelving - multiples depth choices with various retaining lip styles
- Mobile and suspended cabinets - numerous door and drawer configurations
- Horizontal power strip
- Bin box mount




Single Sided Example



Double Sided Example

Agility™ Adaptable  
Table system part numbering guide

AA- 60 30 A 84 R SC



P/N example for 60"W x 30"D x 84"H, adjustable height bench with right-hand Service Hub:  
AA-6030A84R




## Agility Essential Workstations

- Multiple lengths, heights and depths available to meet your specifications
- 11ga steel upright posts
- Knockdown design (assembly required)
- Static load capacity, equally distributed combined system: 3000 lbs


### Available options:

- Pre-wired / pre-plumbed Service Hub system offers numerous power, data and specialty gas services
- 1" thick lab-grade phenolic work surface (Epoxy or stainless also available)
- Upper shelving - multiples depth choices with various retaining lip styles
- Mobile cabinets - numerous door and drawer configurations
- Horizontal power strip
- Bin box mount



## Agility™ Essential Table system part numbering guide

**AE- 48 30 H 84 R**




P/N example for 48"W x 30"D x 84"H, fixed-height (36.75 nom) bench with right-hand Service Hub:  
AE-4830H84R




# AGILITY ESSENTIAL DOUBLE ISLAND CLOSURE PANELS

## Double Island - End View



## Double Island - Front View



## End Filler Panels

- Used to completely close off the space between the leg supports of Agility Essential double island workstations

Model	Table Height
EFP25	29-3/4"
EFP31	35-3/4"


## Rear center Back Filler Panels

- Used to completely close off the back of each side of Agility Essential double island workstations
- Two panes required for double-sided island/benches

Model	Table Height	Table Width
RBP-L4-G2	29-3/4"	47-3/4"
RBP-L5-G2	29-3/4"	59-3/4"
RBP-L6-G2	29-3/4"	71-3/4"
RBP-A4-G2	32"	47-3/4"
RBP-A5-G2	32"	59-3/4"
RBP-A6-G2	32"	71-3/4"
RBP-H4-G2	35-3/4"	47-3/4"
RBP-H5-G2	35-3/4"	59-3/4"
RBP-H6-G2	35-3/4"	71-3/4"

### Agility Service Hub:

Service hubs are an addition that is unique to HLS Agility Workstations. Designed to be used on the Agility Mobile, Agility Adaptable and Agility Essential workstations. This exclusive design distributes power and services vertically for the most efficient use of the workspace.




Service Hub:  
Type 1



Service Hub:  
Type 2



Service Hub:  
Type 3



Service Hub:  
Type 4

# AGILITY ACCESSORIES

## Service Hubs

### Type 1 Service Hub - available options

C1 Detail	C2 Detail	C3 Detail	Power Input
N/A	Single 90° or Dual (180° opposed) needle valve turret*	N/A	N/A


### Type 2 Service Hub - available options

C1 Detail	C2 Detail	C3 Detail	Power Input
Up to four (4) 120v 20A duplexes**	Single 90° or Dual (180° opposed) needle valve turret*	RJ45 CAT 6 ethernet jack (specify quantity)	NEMA L5-20P or L14-20P supply power cord (specify power cord plug type)


\*Specify valve qty (1-3) and service type: Air, Vacuum, Argon, Nitrogen, CDA, CO2

\*\*Specify duplex color (white, red, orange, brown, grey), circuit count (1 or 2) - NOTE that each power circuit requires dedicated supply power cord and plug/receptacle. No more than one power type allowed per channel (ex: std power cannot be combined with emergency power in one channel)

Service Hub - Type 1



Service Hub - Type 2



### Type 3 Service Hub - available options

C1 Detail	C2 Detail	C3 Detail	Power Input
Up to four (4) 120v 20A duplexes*	Single 90° or Dual (180° opposed) needle valve turret**	Up to four (4) 120v 20A duplexes*	NEMA L5-20P or L14-20P supply power cord (specify power cord plug type)


### Type 4 Service Hub - available options

C1 Detail	C2 Detail	C3 Detail	Power Input
Up to four (4) 120v 20A duplexes*	RJ45 CAT 6 ethernet jack (specify quantity)	Up to four (4) 120v 20A duplexes*	NEMA L5-20P or L14-20P supply power cord (specify power cord plug type)


\*Specify duplex color (white, red, orange, brown, grey), circuit count (1 or 2) - NOTE that each power circuit requires dedicated supply power cord and plug/receptacle. No more than one power type allowed per channel (ex: std power cannot be combined with emergency power in one channel)

\*\*Specify valve qty (1-2) and service type: Air, Vacuum, Argon, Nitrogen, CDA, CO2

### Service Hub - Type 3



### Service Hub - Type 4




# AGILITY ACCESSORIES


## Service Hubs

### Needle Valve Turrets

- Each valve includes one (1) pre-plumbed gas hose and hose end per valve
- Specify gas service valve qty (1-3) and service type (Air, Vacuum, Argon, Nitrogen, CDA, CO2)
- Hoses extend out the top of post




90° Valve



180° Valve


### 120V 20A power (specify color)

- Pre-wired 120VAC 20A duplex
- Includes single 12AWG 3C power cable per circuit with L5-20P plug extending out the top of the post
- Specify duplex color (white, red, orange, brown, grey) and qty (1-4)




### 208/220V power (specify type)

- Pre-wired 208VAC L6-20R receptacle
- Includes 12AWG 3C power cable with L6-20P plug extending out the top of the post



### RJ45 data

- Includes prewired RJ45 data jacks (specify quantity)
- Cables extend out the top of post



## M-Series Bench Systems

- HLS M-series bench systems are highly configurable, and offer numerous accessory options designed to meet a variety of project budgets
- This new line of benches has been engineered with features to meet your load and usage demands, including knock-down design for economical transportation
- Power, data, and gas service options can be pre-installed into the frame, resulting in a clean, uncompromised work center
- Services can be upgraded or re-configured in the field with minimal effort

## Bench Configurations:

- M1 - Single-sided, 2x3" 14GA upright post
- M2 - Double-sided, 2x3" 14GA upright post
- M3 - Single-sided, 2x5" 11GA upright post
- M4 - Double-sided, 2x5" 11GA upright post

## Sizing choices


- Widths - 3', 4', 5', 6'
- Depths - 30", 36", 60" (double-sided)

## Included features


- Six-position adjustable work top height (1" increments from 30-3/4 to 36-3/4")
- Knockdown design (assembly required)
- 14 or 11-gauge steel tube construction typical
- Durable powdercoat finish
- 84" tall upright posts
- Accessory port covers allow future upgrades as needed
- Static load capacity, equally distributed combined system: 1200 lbs (SEFA category 4)

## Available options:

- Pre-wired / pre-plumbed power, data and specialty gas integration in vertical posts
- 1" thick lab-grade phenolic work surface (Epoxy or stainless also available)
- Upper shelving - multiple depth choices with various retaining lip styles
- Mobile and suspended cabinets - numerous door and drawer configurations
- Horizontal power strip
- Bin box mount



M1 - (shown w/ accessories)



M2 - (shown w/ accessories)



# M-SERIES WORKSTATIONS

## M-Series Bench system part numbering guide

**M3 - 60 30 H-2P**

**Bench configurations ]**

- M1 - Single-sided, 2x3"
- M2 - Double-sided, 2x3"
- M3 - Single-sided, 2x5"
- M4 - Double-sided, 2x5"

**Work top length ]**

- 36 - 36"
- 48 - 48"
- 60 - 60"
- 72 - 72"

**Work top depth ]**

- 30 - 30" (M1 and M3)
- 36 - 36" (M1 and M3)
- 60 - 60" (M2 and M4)

**Upright post height ]**


- H - 84"

**Optional features ]**


- 1P - 1-post power
- 2P - 2-post power
- 1G2D - 1-gas + 2-data in one post
- 2G2D - 2-gas + 2-data in one post
- 3G2D - 3-gas + 2-data in one post
- 1P1G2D - 1-post power, 1-gas + 2-data other post
- 1P2G2D - 1-post power, 2-gas + 2-data other post
- 1P3G2D - 1-post power, 3-gas + 2-data other pos

\* if no optional features are required, leave part number suffix blank

P/N example for 60"W x 30" M3 frame, with power in both upright posts: M3-6030H-2P



M3 - (shown w/ accessories)



M4 - (shown w/ accessories)

**Application notes-**

M1 & M2: Max. five duplexes per power cord; Max one power cord/circuit per post

M3 & M4: Max. five duplexes per cord; Max two power cords/circuits per post



Part Number	Description	Frame Width	Frame Depth	Top Size (not included)	Frame configuration
<b>M1-6030H-1P1G2D</b>	BENCH FRAME ASY, SINGLE, M1, 60W x 30 X 84H, 2X3, 1-POST PWR, 1-GAS + 2-DATA, NO TOP	59-3/4"	28"	60" x 30"	Single-side
<b>M1-6030H-1P2G2D</b>	BENCH FRAME ASY, SINGLE, M1, 60W x 30 X 84H, 2X3, 1-POST PWR, 2-GAS + 2-DATA, NO TOP	59-3/4"	28"	60" x 30"	Single-side
<b>M1-6030H-1P3G2D</b>	BENCH FRAME ASY, SINGLE, M1, 60W x 30 X 84H, 2X3, 1-POST PWR, 3-GAS + 2-DATA, NO TOP	59-3/4"	28"	60" x 30"	Single-side
<b>M1-7230H-1P1G2D</b>	BENCH FRAME ASY, SINGLE, M1, 72W x 30 X 84H, 2X3, 1-POST PWR, 1-GAS + 2-DATA, NO TOP	71-3/4"	28"	72" x 30"	Single-side
<b>M1-7230H-1P2G2D</b>	BENCH FRAME ASY, SINGLE, M1, 72W x 30 X 84H, 2X3, 1-POST PWR, 2-GAS + 2-DATA, NO TOP	71-3/4"	28"	72" x 30"	Single-side
<b>M1-7230H-1P3G2D</b>	BENCH FRAME ASY, SINGLE, M1, 72W x 30 X 84H, 2X3, 1-POST PWR, 3-GAS + 2-DATA, NO TOP	71-3/4"	28"	72" X 30"	Single-side
<b>M2-6060H-1P1G2D</b>	BENCH FRAME ASY, DOUBLE, M2, 60W x 60 X 84H, 2X3, 1-POST PWR, 1-GAS + 2-DATA, NO TOP	59-3/4"	58"	60" x 60" (2 piece top)	Double-side
<b>M2-6060H-1P2G2D</b>	BENCH FRAME ASY, DOUBLE, M2, 60W x 60 X 84H, 2X3, 1-POST PWR, 2-GAS + 2-DATA, NO TOP	59-3/4"	58"	60" x 60" (2 piece top)	Double-side
<b>M2-6060H-1P3G2D</b>	BENCH FRAME ASY, DOUBLE, M2, 60W x 60 X 84H, 2X3, 1-POST PWR, 3-GAS + 2-DATA, NO TOP	59-3/4"	58"	60" x 60" (2 piece top)	Double-side
<b>M2-7260H-1P1G2D</b>	BENCH FRAME ASY, DOUBLE, M2, 72W x 60 X 84H, 2X3, 1-POST PWR, 1-GAS + 2-DATA, NO TOP	71-3/4"	58"	72" x 60" (2 piece top)	Double-side
<b>M2-7260H-1P2G2D</b>	BENCH FRAME ASY, DOUBLE, M2, 72W x 60 X 84H, 2X3, 1-POST PWR, 2-GAS + 2-DATA, NO TOP	71-3/4"	58"	72" x 60" (2 piece top)	Double-side
<b>M2-7260H-1P3G2D</b>	BENCH FRAME ASY, DOUBLE, M2, 72W x 60 X 84H, 2X3, 1-POST PWR, 3-GAS + 2-DATA, NO TOP	71-3/4"	58"	72" x 60" (2 piece top)	Double-side

NOTE: specify gas service type(s) at time of order.

### Work surfaces; M-Series (Science-grade black phenolic, 1" thick)

Part Number	Description
CTOP-PHENM13630-1.0IN	COUNTERTOP, PHENOLIC, M1 36x30 BLACK, FIN EDGES, 1" THK, EACH
CTOP-PHENM14830-1.0IN	COUNTERTOP, PHENOLIC, M1 48x30 BLACK, FIN EDGES, 1" THK, EACH
CTOP-PHENM16030-1.0IN	COUNTERTOP, PHENOLIC, M1 60x30 BLACK, FIN EDGES, 1" THK, EACH
CTOP-PHENM17230-1.0IN	COUNTERTOP, PHENOLIC, M1 72x30 BLACK, FIN EDGES, 1" THK, EACH
CTOP-PHENM236DS-1.0IN	COUNTERTOP, PHENOLIC, M2 36x30 DS BLACK, FIN EDGES, 1" THK, EA (2 REQ'D)
CTOP-PHENM248DS-1.0IN	COUNTERTOP, PHENOLIC, M2 48x30 DS BLACK, FIN EDGES, 1" THK, EA (2 REQ'D)
CTOP-PHENM260DS-1.0IN	COUNTERTOP, PHENOLIC, M2 60x30 DS BLACK, FIN EDGES, 1" THK, EA (2 REQ'D)
CTOP-PHENM272DS-1.0IN	COUNTERTOP, PHENOLIC, M2 72x30 DS BLACK, FIN EDGES, 1" THK, EA (2 REQ'D)
CTOP-PHENAA3630-1.0IN	COUNTERTOP, PHENOLIC, AA/AE/M3 36x30 BLACK, FIN EDGES, 1" THK, EACH
CTOP-PHENAA4830-1.0IN	COUNTERTOP, PHENOLIC, AA/AE/M3 48x30 BLACK, FIN EDGES, 1" THK, EACH
CTOP-PHENAA6030-1.0IN	COUNTERTOP, PHENOLIC, AA/AE/M3 60x30 BLACK, FIN EDGES, 1" THK, EACH
CTOP-PHENAA7230-1.0IN	COUNTERTOP, PHENOLIC, AA/AE/M3 72x30 BLACK, FIN EDGES, 1" THK, EACH
CTOP-PHENAA36DS-1.0IN	COUNTERTOP, PHENOLIC, AA/AE/M4 36x30 DS BLACK, FIN EDGES, 1" THK, EA (2 REQ'D)
CTOP-PHENAA48DS-1.0IN	COUNTERTOP, PHENOLIC, AA/AE/M4 48x30 DS BLACK, FIN EDGES, 1" THK, EA (2 REQ'D)
CTOP-PHENAA60DS-1.0IN	COUNTERTOP, PHENOLIC, AA/AE/M4 60x30 DS BLACK, FIN EDGES, 1" THK, EA (2 REQ'D)
CTOP-PHENAA72DS-1.0IN	COUNTERTOP, PHENOLIC, AA/AE/M4 72x30 DS BLACK, FIN EDGES, 1" THK, EA (2 REQ'D)


### Suspended Cabinet Mount Kit

(see catalog page 29 for compatible HLS suspended cabinets)


Part Number	Description
K13429-36	MOUNT RAIL KIT, SUSP CAB, 36W, M-SERIES BENCH
K13429-48	MOUNT RAIL KIT, SUSP CAB, 48W, M-SERIES BENCH
K13429-60	MOUNT RAIL KIT, SUSP CAB, 60W, M-SERIES BENCH
K13429-72	MOUNT RAIL KIT, SUSP CAB, 72W, M-SERIES BENCH

### M1/M2 Applications

Style 1A	Style 2A	Style 3A
<p>Pre-wired 120VAC 20A duplexes on single circuit. Includes single 12AWG 3C power cable with L5-20P plug extending out the top of the post. Specify duplex color (white, red, orange, brown, grey) and qty (1-4)</p>	<p>Includes one (1) pre-plumbed gas hose and hose end per valve. Specify gas service valve qty (1-3), service type (Air, Vacuum, Argon, Nitrogen, CDA, CO2), and single (90°) or dual (180°) outlet valve(s). Includes two (2) pre-wired RJ45 ethernet jacks. Hose and data cable extend out the top of post.</p>	<p>Pre-wired 208VAC L6-20R receptacle. Includes 12AWG 3C power cable with L6-20P plug extending out the top of the post.</p>



STYLE 1A



STYLE 2A



STYLE 3A




# M-SERIES


## M3/M4 Power Packs

### M3/M4 Applications

Style 1B	Style 2B	Style 3B	Style 4B	Style 5B
<p>Pre-wired 120VAC 20A duplexes on single circuit. Includes single 12AWG 3C power cable with L5-20P plug extending out the top of the post. Specify duplex color (white, red, orange, brown, grey) and qty (1-4). Two circuits (2 cords) may be available*.</p>	<p>Includes one (1) pre-plumbed gas hose and hose end per valve. Specify gas service valve qty (1-3), service type (Air, Vacuum, Argon, Nitrogen, CDA, CO2), and single (90°) or dual (180°) outlet valve(s). Includes four (4) pre-wired RJ45 ethernet jacks. Hose and data cable extend out the top of post.</p>	<p>Pre-wired 208VAC L6-20R receptacle. Includes 12AWG 3C power cable with L6-20P plug extending out the top of the post.</p>	<p>Pre-wired 120VAC 20A duplexes on single circuit. Includes single 12AWG 3C power cable with L5-20P plug extending out the top of the post. Specify duplex color (white, red, orange, brown, grey) and qty (1-4). Two circuits (2 cords) may be available*. Single 4-port USB power duplex. Integrated 24v 60W power w/ 5.5mm jack for LED undershelf light.</p>	<p>Pre-wired 120VAC 20A duplexes. Includes 12AWG 3C power cable(s) with L5-20P plug(s) extending out the top of the post. Specify duplex color (white, red, orange, brown, grey) and qty (5-8). Two circuits (2 cords) may be available*.</p>




STYLE 1B




STYLE 2B



STYLE 3B



STYLE 4B



STYLE 5B

\* Contact your HLS Sales representative for more details



STYLE 1B




STYLE 2B



STYLE 3B



STYLE 4B





STYLE 5B

# M-SERIES


## Power, Gas and Data Options

### Needle Valve Turrets

- Each valve includes one (1) pre-plumbed gas hose and hose end per valve
- Specify gas service valve qty (1-3) and service type (Air, Vacuum, Argon, Nitrogen, CDA, CO2)
- Hoses extend out the top of post




90° Valve



180° Valve


### 120V 20A power (specify color)

- Pre-wired 120VAC 20A duplexes
- Includes single 12AWG 3C power cable per circuit with L5-20P plug extending out the top of the post
- Specify duplex color (white, red, orange, brown, grey) and qty (1-4)




### 208/220V power (specify type)

- Pre-wired 208VAC L6-20R receptacle
- Includes 12AWG 3C power cable with L6-20P plug extending out the top of the post




### RJ45 data

- Includes prewired RJ45 data jacks (specify quantity)
- Cables extend out the top of post





M1 / M2



M1 / M2




M3 / M4

# TABLES – HEAVY DUTY WELDED

## Fixed Height

- Heavy Duty welded construction, complete with leveling glides
- Tables can be used with a bottom shelf, or mobile cabinet
- Countertop overhangs side edges of table frames by 1/8", 1" overhang at front and rear of frame is typical
- Countertop, brackets, shelves, and cabinets must be ordered separately
- Shipped fully assembled, except for shelves, countertop, or accessories



### Heavy Duty Welded Tables, 29-3/4" High

Model	Frame Width	Frame Depth	Top Size
HDW36A29-G2	35-3/4"	23"	36" x 25"
HDW48A29-G2	47-3/4"	23"	48" x 25"
HDW60A29-G2	59-3/4"	23"	60" x 25"
HDW72A29-G2	71-3/4"	23"	72" x 25"
HDW36B29-G2	35-3/4"	28"	36" x 30"
HDW48B29-G2	47-3/4"	28"	48" X 30"
HDW60B29-G2	59-3/4"	28"	60" X 30"
HDW72B29-G2	71-3/4"	28"	72" X 30"
HDW36C29-G2	35-3/4"	34"	36" X 36"
HDW48C29-G2	47-3/4"	34"	48" X 36"
HDW60C29-G2	59-3/4"	34"	60" X 36"
HDW72C29-G2	71-3/4"	34"	72" X 36"
HDW36D29-G2	35-3/4"	40"	36" X 42"
HDW48D29-G2	47-3/4"	40"	48" x 42"
HDW60D29-G2	59-3/4"	40"	60" x 42"
HDW72D29-G2	71-3/4"	40"	72" x 42"


### Heavy Duty Welded Tables, 35-3/4" High

Model	Frame Width	Frame Depth	Top Size
HDW36A35-G2	35-3/4"	23"	36" x 25"
HDW48A35-G2	47-3/4"	23"	48" x 25"
HDW60A35-G2	59-3/4"	23"	60" x 25"
HDW72A35-G2	71-3/4"	23"	72" x 25"
HDW36B35-G2	35-3/4"	28"	36" x 30"
HDW48B35-G2	47-3/4"	28"	48" X 30"
HDW60B35-G2	59-3/4"	28"	60" X 30"
HDW72B35-G2	71-3/4"	28"	72" X 30"
HDW36C35-G2	35-3/4"	34"	36" X 36"
HDW48C35-G2	47-3/4"	34"	48" X 36"
HDW60C35-G2	59-3/4"	34"	60" X 36"
HDW72C35-G2	71-3/4"	34"	72" X 36"
HDW36D35-G2	35-3/4"	40"	36" X 42"
HDW48D35-G2	47-3/4"	40"	48" x 42"
HDW60D35-G2	59-3/4"	40"	60" x 42"
HDW72D35-G2	71-3/4"	40"	72" x 42"

## Fixed Height w/ Uprights

[hansonlab.com](http://hansonlab.com) • 805-498-3121

- Heavy Duty welded construction, complete with leveling glides
- Tables can be used with a bottom shelf, suspended or mobile cabinet
- Countertop overhangs side edges of table frames by 1/8", 1" overhang at front and rear of frame is typical
- Countertop, brackets, shelves, and cabinets must be ordered separately
- Shipped fully assembled, except for shelves, countertop, or accessories
- Total Height 78" to top of upright posts
- Rear upright posts are 2" wide by 2" deep



Heavy Duty Welded Tables, 29-3/4" High w/ Shelf Upright posts

Model	Frame Width	Frame Depth	Top Size
HDWU36A29-G2	35-3/4"	23"	36" x 25"
HDWU48A29-G2	47-3/4"	23"	48" x 25"
HDWU60A29-G2	59-3/4"	23"	60" x 25"
HDWU72A29-G2	71-3/4"	23"	72" x 25"
HDWU36B29-G2	35-3/4"	28"	36" x 30"
HDWU48B29-G2	47-3/4"	28"	48" X 30"
HDWU60B29-G2	59-3/4"	28"	60" X 30"
HDWU72B29-G2	71-3/4"	28"	72" X 30"
HDWU36C29-G2	35-3/4"	34"	36" X 36"
HDWU48C29-G2	47-3/4"	34"	48" X 36"
HDWU60C29-G2	59-3/4"	34"	60" X 36"
HDWU72C29-G2	71-3/4"	34"	72" X 36"
HDWU36D29-G2	35-3/4"	40"	36" X 42"
HDWU48D29-G2	47-3/4"	40"	48" x 42"
HDWU60D29-G2	59-3/4"	40"	60" x 42"
HDWU72D29-G2	71-3/4"	40"	72" x 42"


Heavy Duty Welded Tables, 35-3/4" High w/ Shelf Upright posts

Model	Frame Width	Frame Depth	Top Size
HDWU36A35-G2	35-3/4"	23"	36" x 25"
HDWU48A35-G2	47-3/4"	23"	48" x 25"
HDWU60A35-G2	59-3/4"	23"	60" x 25"
HDWU72A35-G2	71-3/4"	23"	72" x 25"
HDWU36B35-G2	35-3/4"	28"	36" x 30"
HDWU48B35-G2	47-3/4"	28"	48" X 30"
HDWU60B35-G2	59-3/4"	28"	60" X 30"
HDWU72B35-G2	71-3/4"	28"	72" X 30"
HDWU36C35-G2	35-3/4"	34"	36" X 36"
HDWU48C35-G2	47-3/4"	34"	48" X 36"
HDWU60C35-G2	59-3/4"	34"	60" X 36"
HDWU72C35-G2	71-3/4"	34"	72" X 36"
HDWU36D35-G2	35-3/4"	40"	36" X 42"
HDWU48D35-G2	47-3/4"	40"	48" x 42"
HDWU60D35-G2	59-3/4"	40"	60" x 42"
HDWU72D35-G2	71-3/4"	40"	72" x 42"

# TABLES – HEAVY DUTY WELDED

## Fixed Height w/ Suspended Cabinet Mount

- Heavy Duty welded construction, complete with leveling glides
- Countertops overhang side edges of table frames by 1/8", 1" overhang at front and rear of frame is typical
- Countertop, and cabinets must be ordered separately
- Shipped fully assembled, except for cabinets, countertop, or accessories



### Heavy Duty Welded Tables, 29-3/4" High w/ Shelf Upright Posts & Suspended Cabinet Compatibility

Model	Frame Width	Frame Depth	Top Size
HDW36A29SC-G2	35-3/4"	23"	36" x 25"
HDW48A29SC-G2	47-3/4"	23"	48" x 25"
HDW60A29SC-G2	59-3/4"	23"	60" x 25"
HDW72A29SC-G2	71-3/4"	23"	72" x 25"
HDW36B29SC-G2	35-3/4"	28"	36" x 30"
HDW48B29SC-G2	47-3/4"	28"	48" X 30"
HDW60B29SC-G2	59-3/4"	28"	60" X 30"
HDW72B29SC-G2	71-3/4"	28"	72" X 30"
HDW36C29SC-G2	35-3/4"	34"	36" X 36"
HDW48C29SC-G2	47-3/4"	34"	48" X 36"
HDW60C29SC-G2	59-3/4"	34"	60" X 36"
HDW72C29SC-G2	71-3/4"	34"	72" X 36"
HDW36D29SC-G2	35-3/4"	40"	36" X 42"
HDW48D29SC-G2	47-3/4"	40"	48" x 42"
HDW60D29SC-G2	59-3/4"	40"	60" x 42"
HDW72D29SC-G2	71-3/4"	40"	72" x 42"


### Heavy Duty Welded Tables, 35-3/4" High w/ Shelf Upright Posts & Suspended Cabinet Compatibility

Model	Frame Width	Frame Depth	Top Size
HDW36A35SC-G2	35-3/4"	23"	36" x 25"
HDW48A35SC-G2	47-3/4"	23"	48" x 25"
HDW60A35SC-G2	59-3/4"	23"	60" x 25"
HDW72A35SC-G2	71-3/4"	23"	72" x 25"
HDW36B35SC-G2	35-3/4"	28"	36" x 30"
HDW48B35SC-G2	47-3/4"	28"	48" X 30"
HDW60B35SC-G2	59-3/4"	28"	60" X 30"
HDW72B35SC-G2	71-3/4"	28"	72" X 30"
HDW36C35SC-G2	35-3/4"	34"	36" X 36"
HDW48C35SC-G2	47-3/4"	34"	48" X 36"
HDW60C35SC-G2	59-3/4"	34"	60" X 36"
HDW72C35SC-G2	71-3/4"	34"	72" X 36"
HDW36D35SC-G2	35-3/4"	40"	36" X 42"
HDW48D35SC-G2	47-3/4"	40"	48" x 42"
HDW60D35SC-G2	59-3/4"	40"	60" x 42"
HDW72D35SC-G2	71-3/4"	40"	72" x 42"

## Fixed Height w/ Upright Posts & Suspended Cabinet Mount

[hansonlab.com](http://hansonlab.com) • 805-498-3121

- Heavy Duty welded construction, complete with leveling glides
- Countertops overhang side edges of table frames by 1/8", 1" overhang at front and rear of frame is typical
- Countertop, and cabinets must be ordered separately
- Shipped fully assembled, except for cabinets, countertop, or accessories
- Total Height 78" to top of upright posts
- Rear upright posts are 2" wide by 2" deep



### Heavy Duty Welded Tables, 29-3/4" High w/ Shelf Upright Posts & Suspended Cabinet Compatibility

Model	Frame Width	Frame Depth	Top Size
HDWU36A29SC-G2	35-3/4"	23"	36" x 25"
HDWU48A29SC-G2	47-3/4"	23"	48" x 25"
HDWU60A29SC-G2	59-3/4"	23"	60" x 25"
HDWU72A29SC-G2	71-3/4"	23"	72" x 25"
HDWU36B29SC-G2	35-3/4"	28"	36" x 30"
HDWU48B29SC-G2	47-3/4"	28"	48" X 30"
HDWU60B29SC-G2	59-3/4"	28"	60" X 30"
HDWU72B29SC-G2	71-3/4"	28"	72" X 30"
HDWU36C29SC-G2	35-3/4"	34"	36" X 36"
HDWU48C29SC-G2	47-3/4"	34"	48" X 36"
HDWU60C29SC-G2	59-3/4"	34"	60" X 36"
HDWU72C29SC-G2	71-3/4"	34"	72" X 36"
HDWU36D29SC-G2	35-3/4"	40"	36" X 42"
HDWU48D29SC-G2	47-3/4"	40"	48" x 42"
HDWU60D29SC-G2	59-3/4"	40"	60" x 42"
HDWU72D29SC-G2	71-3/4"	40"	72" x 42"


### Heavy Duty Welded Tables, 35-3/4" High w/ Shelf Upright Posts & Suspended Cabinet Compatibility

Model	Frame Width	Frame Depth	Top Size
HDWU36A35SC-G2	35-3/4"	23"	36" x 25"
HDWU48A35SC-G2	47-3/4"	23"	48" x 25"
HDWU60A35SC-G2	59-3/4"	23"	60" x 25"
HDWU72A35SC-G2	71-3/4"	23"	72" x 25"
HDWU36B35SC-G2	35-3/4"	28"	36" x 30"
HDWU48B35SC-G2	47-3/4"	28"	48" X 30"
HDWU60B35SC-G2	59-3/4"	28"	60" X 30"
HDWU72B35SC-G2	71-3/4"	28"	72" X 30"
HDWU36C35SC-G2	35-3/4"	34"	36" X 36"
HDWU48C35SC-G2	47-3/4"	34"	48" X 36"
HDWU60C35SC-G2	59-3/4"	34"	60" X 36"
HDWU72C35SC-G2	71-3/4"	34"	72" X 36"
HDWU36D35SC-G2	35-3/4"	40"	36" X 42"
HDWU48D35SC-G2	47-3/4"	40"	48" x 42"
HDWU60D35SC-G2	59-3/4"	40"	60" x 42"
HDWU72D35SC-G2	71-3/4"	40"	72" x 42"

# TABLES – HEAVY DUTY WELDED

## Fixed Height w/ Casters

- Heavy Duty welded construction, complete with caster kits
- Tables can be used with a bottom shelf, or mobile cabinet
- Countertops overhang side edges of table frames by 1/8", 1" overhang at front and rear of frame is typical
- Countertop, brackets, shelves, and cabinets must be ordered separately
- Shipped fully assembled, except for shelves, countertop, or accessories



Heavy Duty Welded Tables w/ Casters,  
29-3/4" High

Model	Frame Width	Frame Depth	Top Size
HDW36A29C-G2	35-3/4"	23"	36" x 25"
HDW48A29C-G2	47-3/4"	23"	48" x 25"
HDW60A29C-G2	59-3/4"	23"	60" x 25"
HDW72A29C-G2	71-3/4"	23"	72" x 25"
HDW36B29C-G2	35-3/4"	28"	36" x 30"
HDW48B29C-G2	47-3/4"	28"	48" x 30"
HDW60B29C-G2	59-3/4"	28"	60" x 30"
HDW72B29C-G2	71-3/4"	28"	72" x 30"
HDW36C29C-G2	35-3/4"	34"	36" x 36"
HDW48C29C-G2	47-3/4"	34"	48" x 36"
HDW60C29C-G2	59-3/4"	34"	60" x 36"
HDW72C29C-G2	71-3/4"	34"	72" x 36"
HDW36D29C-G2	35-3/4"	40"	36" x 42"
HDW48D29C-G2	47-3/4"	40"	48" x 42"
HDW60D29C-G2	59-3/4"	40"	60" x 42"
HDW72D29C-G2	71-3/4"	40"	72" x 42"


Heavy Duty Welded Tables w/  
Casters, 35-3/4" High

Model	Frame Width	Frame Depth	Top Size
HDW36A35C-G2	35-3/4"	23"	36" x 25"
HDW48A35C-G2	47-3/4"	23"	48" x 25"
HDW60A35C-G2	59-3/4"	23"	60" x 25"
HDW72A35C-G2	71-3/4"	23"	72" x 25"
HDW36B35C-G2	35-3/4"	28"	36" x 30"
HDW48B35C-G2	47-3/4"	28"	48" x 30"
HDW60B35C-G2	59-3/4"	28"	60" x 30"
HDW72B35C-G2	71-3/4"	28"	72" x 30"
HDW36C35C-G2	35-3/4"	34"	36" x 36"
HDW48C35C-G2	47-3/4"	34"	48" x 36"
HDW60C35C-G2	59-3/4"	34"	60" x 36"
HDW72C35C-G2	71-3/4"	34"	72" x 36"
HDW36D35C-G2	35-3/4"	40"	36" x 42"
HDW48D35C-G2	47-3/4"	40"	48" x 42"
HDW60D35C-G2	59-3/4"	40"	60" x 42"
HDW72D35C-G2	71-3/4"	40"	72" x 42"

## Fixed Height w/ Upright Posts & Casters

[hansonlab.com](http://hansonlab.com) • 805-498-3121

- Heavy Duty welded construction, complete with caster kits
- Tables can be used with a bottom shelf, or mobile cabinet
- Countertops overhang side edges of table frames by 1/8", 1" overhang at front and rear of frame is typical
- Countertop, brackets, shelves, and cabinets must be ordered separately
- Shipped fully assembled, except for shelves, countertop, or accessories
- Total Height 78" to top of upright posts
- Rear upright posts are 2" wide by 2" deep



### Heavy Duty Welded Tables w/ Upright Posts & Casters, 29-3/4" High

Model	Frame Width	Frame Depth	Top Size
HDWU36A29C-G2	35-3/4"	23"	36" x 25"
HDWU48A29C-G2	47-3/4"	23"	48" x 25"
HDWU60A29C-G2	59-3/4"	23"	60" x 25"
HDWU72A29C-G2	71-3/4"	23"	72" x 25"
HDWU36B29C-G2	35-3/4"	28"	36" x 30"
HDWU48B29C-G2	47-3/4"	28"	48" X 30"
HDWU60B29C-G2	59-3/4"	28"	60" X 30"
HDWU72B29C-G2	71-3/4"	28"	72" X 30"
HDWU36C29C-G2	35-3/4"	34"	36" X 36"
HDWU48C29C-G2	47-3/4"	34"	48" X 36"
HDWU60C29C-G2	59-3/4"	34"	60" X 36"
HDWU72C29C-G2	71-3/4"	34"	72" X 36"
HDWU36D29C-G2	35-3/4"	40"	36" X 42"
HDWU48D29C-G2	47-3/4"	40"	48" x 42"
HDWU60D29C-G2	59-3/4"	40"	60" x 42"
HDWU72D29C-G2	71-3/4"	40"	72" x 42"


### Heavy Duty Welded Tables w/ Upright Posts & Casters, 35-3/4" High

Model	Frame Width	Frame Depth	Top Size
HDWU36A35C-G2	35-3/4"	23"	36" x 25"
HDWU48A35C-G2	47-3/4"	23"	48" x 25"
HDWU60A35C-G2	59-3/4"	23"	60" x 25"
HDWU72A35C-G2	71-3/4"	23"	72" x 25"
HDWU36B35C-G2	35-3/4"	28"	36" x 30"
HDWU48B35C-G2	47-3/4"	28"	48" X 30"
HDWU60B35C-G2	59-3/4"	28"	60" X 30"
HDWU72B35C-G2	71-3/4"	28"	72" X 30"
HDWU36C35C-G2	35-3/4"	34"	36" X 36"
HDWU48C35C-G2	47-3/4"	34"	48" X 36"
HDWU60C35C-G2	59-3/4"	34"	60" X 36"
HDWU72C35C-G2	71-3/4"	34"	72" X 36"
HDWU36D35C-G2	35-3/4"	40"	36" X 42"
HDWU48D35C-G2	47-3/4"	40"	48" x 42"
HDWU60D35C-G2	59-3/4"	40"	60" x 42"
HDWU72D35C-G2	71-3/4"	40"	72" x 42"

# TABLES – HEAVY DUTY WELDED


## Adjustable Height

- Heavy Duty welded construction, complete with leveling glides
- Tables can be used with a bottom shelf, or mobile cabinet
- Countertops overhang side edges of table frames by 1/8", 1" overhang at front and rear of frame is typical
- Countertop, brackets, shelves, and cabinets must be ordered separately
- Shipped fully assembled, except for shelves, countertop, or accessories
- Total Height 78" to top of upright posts
- Rear upright posts are 2" wide by 2" deep
- Frame height adjustment is from 29-3/4" - 35-3/4" w/ out countertop



### Heavy Duty Welded Tables Adjustable Height

Model	Frame Width	Frame Depth	Top Size
HDW36AADJ-G2	35-3/4"	23"	36" x 25"
HDW48AADJ-G2	47-3/4"	23"	48" x 25"
HDW60AADJ-G2	59-3/4"	23"	60" x 25"
HDW72AADJ-G2	71-3/4"	23"	72" x 25"
HDW36BADJ-G2	35-3/4"	28"	36" x 30"
HDW48BADJ-G2	47-3/4"	28"	48" X 30"
HDW60BADJ-G2	59-3/4"	28"	60" X 30"
HDW72BADJ-G2	71-3/4"	28"	72" X 30"
HDW36CADJ-G2	35-3/4"	34"	36" X 36"
HDW48CADJ-G2	47-3/4"	34"	48" X 36"
HDW60CADJ-G2	59-3/4"	34"	60" X 36"
HDW72CADJ-G2	71-3/4"	34"	72" X 36"
HDW36DADJ-G2	35-3/4"	40"	36" X 42"
HDW48DADJ-G2	47-3/4"	40"	48" x 42"
HDW60DADJ-G2	59-3/4"	40"	60" x 42"
HDW72DADJ-G2	71-3/4"	40"	72" x 42"




### Heavy Duty Welded Tables, Adjustable Height w/ Shelf Upright Posts

Model	Frame Width	Frame Depth	Top Size
HDWU36AADJ-G2	35-3/4"	23"	36" x 25"
HDWU48AADJ-G2	47-3/4"	23"	48" x 25"
HDWU60AADJ-G2	59-3/4"	23"	60" x 25"
HDWU72AADJ-G2	71-3/4"	23"	72" x 25"
HDWU36BADJ-G2	35-3/4"	28"	36" x 30"
HDWU48BADJ-G2	47-3/4"	28"	48" X 30"
HDWU60BADJ-G2	59-3/4"	28"	60" X 30"
HDWU72BADJ-G2	71-3/4"	28"	72" X 30"
HDWU36CADJ-G2	35-3/4"	34"	36" X 36"
HDWU48CADJ-G2	47-3/4"	34"	48" X 36"
HDWU60CADJ-G2	59-3/4"	34"	60" X 36"
HDWU72CADJ-G2	71-3/4"	34"	72" X 36"
HDWU36DADJ-G2	35-3/4"	40"	36" X 42"
HDWU48DADJ-G2	47-3/4"	40"	48" x 42"
HDWU60DADJ-G2	59-3/4"	40"	60" x 42"
HDWU72DADJ-G2	71-3/4"	40"	72" x 42"

## Adjustable Height w/ Casters

[hansonlab.com](http://hansonlab.com) • 805-498-3121


- Heavy Duty welded construction, complete with casters
- Tables can be used with a bottom shelf, or mobile cabinet
- Countertops overhang side edges of table frames by 1/8", 1" overhang at front and rear of frame is typical
- Countertop, brackets, shelves, and cabinets must be ordered separately
- Shipped fully assembled, except for shelves, countertop, or accessories
- Total Height 78" to top of upright posts
- Rear upright posts are 2" wide by 2" deep
- Frame height adjustment is from 29-3/4" - 35-3/4" w/ out countertop



### Heavy Duty Welded Tables, Adjustable Height w/ Casters

Model	Frame Width	Frame Depth	Top Size
HDW36AADJC-G2	35-3/4"	23"	36" x 25"
HDW48AADJC-G2	47-3/4"	23"	48" x 25"
HDW60AADJC-G2	59-3/4"	23"	60" x 25"
HDW72AADJC-G2	71-3/4"	23"	72" x 25"
HDW36BADJC-G2	35-3/4"	28"	36" x 30"
HDW48BADJC-G2	47-3/4"	28"	48" x 30"
HDW60BADJC-G2	59-3/4"	28"	60" x 30"
HDW72BADJC-G2	71-3/4"	28"	72" x 30"
HDW36CADJC-G2	35-3/4"	34"	36" x 36"
HDW48CADJC-G2	47-3/4"	34"	48" x 36"
HDW60CADJC-G2	59-3/4"	34"	60" x 36"
HDW72CADJC-G2	71-3/4"	34"	72" x 36"
HDW36DADJC-G2	35-3/4"	40"	36" x 42"
HDW48DADJC-G2	47-3/4"	40"	48" x 42"
HDW60DADJC-G2	59-3/4"	40"	60" x 42"
HDW72DADJC-G2	71-3/4"	40"	72" x 42"

### Heavy Duty Welded Tables, Adjustable Height w/ Shelf Upright Posts & Casters




Model	Frame Width	Frame Depth	Top Size
HDWU36AADJC-G2	35-3/4"	23"	36" x 25"
HDWU48AADJC-G2	47-3/4"	23"	48" x 25"
HDWU60AADJC-G2	59-3/4"	23"	60" x 25"
HDWU72AADJC-G2	71-3/4"	23"	72" x 25"
HDWU36BADJC-G2	35-3/4"	28"	36" x 30"
HDWU48BADJC-G2	47-3/4"	28"	48" x 30"
HDWU60BADJC-G2	59-3/4"	28"	60" x 30"
HDWU72BADJC-G2	71-3/4"	28"	72" x 30"
HDWU36CADJC-G2	35-3/4"	34"	36" x 36"
HDWU48CADJC-G2	47-3/4"	34"	48" x 36"
HDWU60CADJC-G2	59-3/4"	34"	60" x 36"
HDWU72CADJC-G2	71-3/4"	34"	72" x 36"
HDWU36DADJC-G2	35-3/4"	40"	36" x 42"
HDWU48DADJC-G2	47-3/4"	40"	48" x 42"
HDWU60DADJC-G2	59-3/4"	40"	60" x 42"
HDWU72DADJC-G2	71-3/4"	40"	72" x 42"

# TABLES – HEAVY DUTY WELDED


## Adjustable Height w/ Suspended Cabinet Mount & Casters

- Heavy Duty welded construction, complete with casters
- Countertops overhang side edges of table frames by 1/8", 1" overhang at front and rear of frame is typical
- Countertop, brackets, shelves, and cabinets must be ordered separately
- Shipped fully assembled, except for shelves, countertop, or accessories
- Total Height 78" to top of upright posts
- Rear upright posts are 2" wide by 2" deep
- Frame height adjustment is from 29-3/4" - 35-3/4" w/ out countertop



### Heavy Duty Welded Tables, Adjustable Height w/ Suspended Cabinet Compatibility & Casters

Model	Frame Width	Frame Depth	Top Size
HDW36ASCADJC-G2	35-3/4"	23"	36" x 25"
HDW48ASCADJC-G2	47-3/4"	23"	48" x 25"
HDW60ASCADJC-G2	59-3/4"	23"	60" x 25"
HDW72ASCADJC-G2	71-3/4"	23"	72" x 25"
HDW36BSCADJC-G2	35-3/4"	28"	36" x 30"
HDW48BSCADJC-G2	47-3/4"	28"	48" X 30"
HDW60BSCADJC-G2	59-3/4"	28"	60" X 30"
HDW72BSCADJC-G2	71-3/4"	28"	72" X 30"
HDW36CSCADJC-G2	35-3/4"	34"	36" X 36"
HDW48CSCADJC-G2	47-3/4"	34"	48" X 36"
HDW60CSCADJC-G2	59-3/4"	34"	60" X 36"
HDW72CSCADJC-G2	71-3/4"	34"	72" X 36"
HDW36DSCADJC-G2	35-3/4"	40"	36" X 42"
HDW48DSCADJC-G2	47-3/4"	40"	48" x 42"
HDW60DSCADJC-G2	59-3/4"	40"	60" x 42"
HDW72DSCADJC-G2	71-3/4"	40"	72" x 42"




### Heavy Duty Welded Tables, Adjustable Height w/ Shelf Upright Posts, Suspended Cabinet Compatibility & Casters

Model	Frame Width	Frame Depth	Top Size
HDWU36ASCADJC-G2	35-3/4"	23"	36" x 25"
HDWU48ASCADJC-G2	47-3/4"	23"	48" x 25"
HDWU60ASCADJC-G2	59-3/4"	23"	60" x 25"
HDWU72ASCADJC-G2	71-3/4"	23"	72" x 25"
HDWU36BSCADJC-G2	35-3/4"	28"	36" x 30"
HDWU48BSCADJC-G2	47-3/4"	28"	48" X 30"
HDWU60BSCADJC-G2	59-3/4"	28"	60" X 30"
HDWU72BSCADJC-G2	71-3/4"	28"	72" X 30"
HDWU36CSCADJC-G2	35-3/4"	34"	36" X 36"
HDWU48CSCADJC-G2	47-3/4"	34"	48" X 36"
HDWU60CSCADJC-G2	59-3/4"	34"	60" X 36"
HDWU72CSCADJC-G2	71-3/4"	34"	72" X 36"
HDWU36DSCADJC-G2	35-3/4"	40"	36" X 42"
HDWU48DSCADJC-G2	47-3/4"	40"	48" x 42"
HDWU60DSCADJC-G2	59-3/4"	40"	60" x 42"
HDWU72DSCADJC-G2	71-3/4"	40"	72" x 42"

- Tables can be used with upright posts, upper shelves, lower shelf, power strip kit, single drawer or mobile cabinet
- Countertop, brackets, shelves, and cabinets must be ordered separately
- Typical countertop overhang on Utility Table frame is 1/8" at sides, 1" at front, flush at rear (for optional upright kit installation)
- Assembly required

### Adjustable Height Frames


Frame height adjustment is from 29-3/4" - 35-3/4" w/out countertop



Model	Frame Width	Frame Depth	Top Size
UTA36B	35-3/4"	29"	36" x 30"
UTA48B	47-3/4"	29"	48" x 30"
UTA60B	59-3/4"	29"	60" x 30"
UTA72B	71-3/4"	29"	72" x 30"
UTA36C	35-3/4"	35"	36" x 36"
UTA48C	47-3/4"	35"	48" X 36"
UTA60C	59-3/4"	35"	60" X 36"
UTA72C	71-3/4"	35"	72" X 36"

### Fixed Height Frames


Frame height is 35-3/4" w/out countertop



Model	Frame Width	Frame Depth	Top Size
UTFH36B	35-3/4"	29"	36" x 30"
UTFH48B	47-3/4"	29"	48" x 30"
UTFH60B	59-3/4"	29"	60" x 30"
UTFH72B	71-3/4"	29"	72" x 30"
UTFH36C	35-3/4"	35"	36" x 36"
UTFH48C	47-3/4"	35"	48" X 36"
UTFH60C	59-3/4"	35"	60" X 36"
UTFH72C	71-3/4"	35"	72" X 36"

### Fixed Height Frames


Frame height is 29-3/4" w/out countertop



Model	Frame Width	Frame Depth	Top Size
UTFL36B	35-3/4"	29"	36" x 30"
UTFL48B	47-3/4"	29"	48" x 30"
UTFL60B	59-3/4"	29"	60" x 30"
UTFL72B	71-3/4"	29"	72" x 30"
UTFL36C	35-3/4"	35"	36" x 36"
UTFL48C	47-3/4"	35"	48" X 36"
UTFL60C	59-3/4"	35"	60" X 36"
UTFL72C	71-3/4"	35"	72" X 36"

# UTILITY TABLES


## Accessories



### Upright Post Kits

Includes 2" X 2" vertical posts, horizontal cross members, and hardware. Table, shelves, and shelf brackets not included


Model	Frame Width	Frame Height
K13409-36	35-3/4"	52-1/2"
K13409-48	47-3/4"	52-1/2"
K13409-60	59-3/4"	52-1/2"
K13409-72	71-3/4"	52-1/2"



### Lower Shelving Kits

Under counter use includes shelf, & mounting brackets. Table frame, & work surface sold separately.


Model	Length	Depth
K7003-36	35-1/2"	24-3/4"
K7003-48	47-1/2"	24-3/4"
K7003-60	59-3/4"	24-3/4"
K7003-72	71-3/4"	24-3/4"



### Rear Apron

Suggested for applications with dynamic equipment


Model	Length	Height
C13406-36	36"	6-1/2"
C13406-48	48"	6-1/2"
C13406-60	60"	6-1/2"
C13406-72	72"	6-1/2"



### Frame Depth Extension Apron

Frame depth extension aprons accommodates 24" or 36" deep work surface (sold separately)

Model	Apron Width	Finished Frame Depth
K13408-24-1	18"	22"
K13408-24	19"	23"
K13408-36-1	30"	34"
K13408-36	31"	35"




### Single Drawer

Includes single 24" x 28" drawer housing. Table frame and work surface sold separately.

Model	Width	Depth
		28"

## Work Surface


HDW and Agility benches w/ out upright posts and all Utility table applications. Chemical resistant black phenolic resin (36" deep sizes can be custom ordered)



Model	Length	Depth	Thickness
CTOP-PHEN3630-1.0IN	36"	30"	1"
CTOP-PHEN4830-1.0IN	48"	30"	1"
CTOP-PHEN6030-1.0IN	60"	30"	1"
CTOP-PHEN7230-1.0IN	72"	30"	1"

## Caster Kit



1/2"-13 Stem swivel casters with brake, use with HDW & Utility Table for light duty applications



Model	Size
H260	Caster, 3"Ø Lt Med Gray 4.25" OAH
H261	Caster, 3.5"Ø Lt Med Gray 4.75" OAH
H262	Caster, 4"Ø Lt Med Gray, 5.25" OAH

## Lower Shelving for Agility Essential & HDW


Painted steel shelf - attaches to lower horizontal support. Table cross brace is suggested for heavy duty applications



Model	Shelf Width	Shelf Depth
TS36A-G2	35-1/2"	18-3/4"
TS36B-G2	35-1/2"	23-1/2"
TS48A-G2	47-1/2"	18-3/4"
TS48B-G2	47-1/2"	23-1/2"
TS60A-G2	59-1/2"	18-3/4"
TS60B-G2	59-1/2"	23-1/2"
TS72A-G2	71-1/2"	18-3/4"
TS72B-G2	71-1/2"	23-1/2"

## Table Cross Brace


Cross brace for added support of the bottom shelf



Model	Brace Width	Work Surface Width
CB36-G2	31-3/4"	36"
CB48-G2	43-3/4"	48"
CB60-G2	55-3/4"	60"
CB72-G2	67-3/4"	72"

# TABLE ACCESSORIES


## Power Strips, Bin Box Mounts, Upper Shelves



### Power Strip Kits (15A)

Upright Post Kits Sold Separately


Model	Mount Width	Strip Length	Cord Length	Outlet Qty.
K13411-36	34.65"	24"	6'	8
K13411-48	46.65"	36"	15'	12
K13411-60	58.65"	48"	15'	8
K13411-72	70.65"	60"	15'	20



### Bin Box Mounts

Bin box not included

Model	Length
A13412-36	36"
A13412-48	48"
A13412-60	60"
A13412-72	72"



### Upper Shelf Kits\*

Includes 12" deep painted metal shelf, and brackets.

Requires upright post kit, and table sold separately.

Model	Width	Depth
K13070-12	31-1/2"	12"
K13032-12	43-1/2"	12"
K13033-12	55-1/2"	12"
K13034-12	67-1/2"	12"

\*see page 103 for full range of HLS shelving sizes & styles