



M1 Single-Sided




M2 Double-Sided


The highly configurable **M-Series** workstations are the ideal solution for the modern lab. These table and bench systems offer numerous accessory options designed to meet a variety of project budgets.

Model	Configuration	Width	Depth	Height
M1	Single-Sided	48", 60" or 72"	30" or 36"	Adjustable up to 1" increments (30-36")
M2	Double-Sided	48", 60" or 72"	30", 36" or 60"	

Power, Data, And Gas Service Options Can Be Pre-installed



Up to four (4) pre-plumbed needle valve service fixtures









Up to four (4) pre-wired duplex receptacles

Seven Durable Powder Coat Finish Options



# HLS M1 & M2 WORKSTATIONS

- **FRAME CONFIGURATIONS**
  - M1: SINGLE-SIDED, 30 OR 36" DEEP
  - M2: DOUBLE-SIDED, 60" DEEP
- **14GA STEEL 2 x 3" REAR POST, 2 x 2" FRONT LEG CONSTRUCTION**
- **SIX-POSITION ADJUSTABLE WORKSURFACE HEIGHT (1" INCREMENTS)**
- **STATIC LOAD CAPACITY, EQUALLY DISTRIBUTED**
  - COMBINED SYSTEM : 1200 LBS
  - LOWER SHELF : 200 LBS
  - UPPER SHELVES : 150 LBS
- **AVAILABLE OPTIONS:**
  - UPPER SHELF KIT (SHOWN) IN VARIOUS SIZES / STYLES
  - LOWER SHELF KIT (SHOWN)
  - POWER, DATA AND GAS SERVICE INTEGRATION INTO UPRIGHT POSTS (SHOWN)
  - CHEM RESISTANT PHENOLIC WORK SURFACE (SHOWN)
  - SIDE APRON OPTION TO ACCOMMODATE ALTERNATE DEPTH WORK SURFACES
  - SUSPENDED CABINET (SHOWN) IN VARIOUS SIZES / STYLES
  - MOBILE CABINET IN VARIOUS SIZES / STYLES
  - HORIZONTAL POWER STRIP KIT (SHOWN)
  - BIN BOX MOUNT BRACKET (SHOWN)
- SEE PRICE LIST FOR FRAME AND ACCESSORY PART NUMBERS
- ASSEMBLY REQUIRED



SIZE MATRIX		
Nominal (width x depth)	Frame width	Worksurface width
36 x 30, 36 or 60*	35-3/4"	36"
48 x 30, 36 or 60*	47-3/4"	48"
60 x 30, 36 or 60*	59-3/4"	60"
72 x 30, 36 or 60*	71-3/4"	72"

\*58" or 60" depth indicates double-sided M2 bench

DRAWN	MWR	4/24/2023	<b>HLS™</b> HANSON LAB SOLUTIONS	
CHECKED			747 Calle Plano, Camarillo CA 93012   805-498-3121	
QA			TITLE	BENCH ASSY, HLS M1/M2 SERIES 2x3, ADJ HT
APPROVED			FINISH	POWDERCOAT PER JOB REQUIREMENT
UNLESS OTHERWISE NOTED, GEN. TOLERANCE +/- .015"; ALL DIMENSIONS ARE IN INCHES	JOB #		SIZE	PART NO
	QUOTE #		D	SEE CATALOG
			SCALE	.13 MATL
				REV A
				SHEET 1 OF 1