



M3 Single-Sided




M4 Double-Sided

The highly configurable **M-Series** workstations are the ideal solution for the modern lab. These table and bench systems offer numerous accessory options designed to meet a variety of project budgets.

| Model | Configuration | Width | Depth | Height |
|-------|---------------|-----------------|------------|---|
| M3 | Single-Sided | 48", 60" or 72" | 30" or 36" | Adjustable up to 1" increments (30-36") |
| M4 | Double-Sided | 48", 60" or 72" | 60" or 72" | |

Power, Data, And Gas Service Options Can Be Pre-installed








Seven Durable Powder Coat Finish Options



HLS M3 & M4 WORKSTATIONS

- **FRAME CONFIGURATIONS**
 - M3: SINGLE-SIDED, 30 OR 36" DEEP
 - M4: DOUBLE-SIDED, 60" DEEP
- **11 GA STEEL 2 x 5" REAR POST, 14 GA 2 x 2" FRONT LEG CONSTRUCTION**
- **SIX-POSITION ADJUSTABLE WORKSURFACE HEIGHT (1" INCREMENTS)**
- **STATIC LOAD CAPACITY, EQUALLY DISTRIBUTED**
 - COMBINED SYSTEM : 1200 LBS
 - LOWER SHELF : 200 LBS
 - UPPER SHELVES : 150 LBS
- **AVAILABLE OPTIONS:**
 - UPPER SHELF KIT (SHOWN) IN VARIOUS SIZES / STYLES
 - LOWER SHELF KIT
 - POWER, DATA AND GAS SERVICE INTEGRATION INTO UPRIGHT POSTS (SHOWN)
 - CHEM RESISTANT PHENOLIC WORK SURFACE (SHOWN)
 - SIDE APRON OPTION TO ACCOMMODATE ALTERNATE DEPTH WORK SURFACES
 - SUSPENDED CABINET (SHOWN) IN VARIOUS SIZES / STYLES
 - MOBILE CABINET IN VARIOUS SIZES / STYLES
 - HORIZONTAL POWER STRIP KIT (SHOWN)
 - BIN BOX MOUNT BRACKET (SHOWN)
- SEE PRICE LIST FOR FRAME AND ACCESSORY PART NUMBERS
- ASSEMBLY REQUIRED



| SIZE MATRIX | | |
|-------------------------|-------------|-------------------|
| Nominal (width x depth) | Frame width | Worksurface width |
| 36 x 30, 36 or 60* | 35-3/4" | 36" |
| 48 x 30, 36 or 60* | 47-3/4" | 48" |
| 60 x 30, 36 or 60* | 59-3/4" | 60" |
| 72 x 30, 36 or 60* | 71-3/4" | 72" |

*58" or 60" depth indicates double-sided M4 bench

| | | | |
|---|------------------|--|------------------------|
| DRAWN MWR | 4/24/2023 | HLS™ HANSON LAB SOLUTIONS 747 Calle Plano, Camarillo CA 93012 805-498-3121 | |
| CHECKED | | TITLE BENCH ASSY, HLS M3 / M4 SERIES 2x5, ADJ HT | |
| QA | | FINISH POWDERCOAT PER JOB REQUIREMENT | |
| APPROVED | | SIZE D | PART NO SEE CATALOG |
| UNLESS OTHERWISE NOTED, GEN. TOLERANCE +/- .015"; ALL DIMENSIONS ARE IN INCHES | JOB # QUOTE # | SCALE .13 | REV A |
| | | MATL | SHEET 1 OF 1 |