

# STEEL CASEWORK

Steel casework is the backbone of your lab.  
Nobody does it better than Hanson.


## Steel Casework Index

	Page
<b>STEEL CABINET DETAILS</b> .....	11
<b>STANDING HEIGHT CABINETS</b> .....	13
• Base Cabinets .....	13
• Corner Base Cabinets .....	17
• Sink Cabinets .....	18
<b>SITTING HEIGHT CABINETS</b> .....	19
• Base Cabinets .....	19
• Corner Base Cabinets .....	23
• Sink Cabinets .....	24
<b>ADA HEIGHT CABINETS</b> .....	25
• Base Cabinets .....	25
• Corner Base Cabinets .....	27
• Sink Cabinets .....	28
<b>SUSPENDED CABINETS</b> .....	29
<b>MOBILE CABINETS</b> .....	30
<b>CHEMICAL STORAGE CABINETS</b> .....	32
• Acid/Corrosive Cabinets .....	32
• Flammable Cabinets.....	35
<b>SPECIALITY CABINETS</b> .....	38
• Equipment Cabinets .....	38
• Vacuum Cabinets.....	39


	Page
<b>CASEWORK ACCESSORIES, FILLERS, APRONS</b> .....	42
• Cabinet Accessories .....	42
• Fillers and Scribes .....	43
• Knee Space Panels .....	44
• Rear Enclosure Panels .....	46
• Drawer Frames, Apron Rails .....	47
• Leg Assemblies and Countertop Support .....	48
• Umbilical Pipe Chases.....	49
<b>WALL CABINETS</b> .....	50
• Wall Cabinet Details .....	50
• 24" Wall Cabinet .....	51
• 30" Wall Cabinet .....	53
• 36" Wall Cabinet .....	55
• 48" Wall Cabinet .....	57
• Wall Cabinet Accessories .....	59
<b>TALL CABINETS</b> .....	61
• Tall Cabinet Details.....	61
• 84" Tall Cabinets.....	62
• Tall Cabinet Accessories.....	67

*\*Laminate and wood may also be available. Contact your HLS Sales representative for more details*





## HLS Features and Benefits



- Ⓐ Cabinet body design accepts multiple arrangements of doors and drawers to permit field changes at any time.
- Ⓑ One piece drawer bodies with coved horizontal corners for easy cleaning.
- Ⓒ 18 ga. double walled drawer and door fronts, with sound deadening full channeled edges.
- Ⓓ Full extension, ball bearing, 100 lb. capacity drawer slides.
- Ⓔ 14 ga. welded gusset plates with hex head leveling bolts at all corners.
- Ⓕ 18 ga. fully removable back panels included in door and open cabinets for access to service piping behind cabinets.
- Ⓖ 20 ga. cabinet bodies with flush seamless welded face. Channel - shaped intermediate uprights prevent cabinet racking and deliver a precise fit of doors and drawers.
- Ⓗ 20 ga. adjustable shelves are included in open and cupboard cabinets.
- Ⓘ Stainless steel 5-knuckle hinges.
- Ⓝ HLS cabinets have electrostatically applied baked on powder coated epoxy paint finish for superior performance.

**THE SCIENCE OF  
LABORATORY  
CONSTRUCTION...**



# STEEL CABINET DETAILS

## General Dimensions

### Hanson Base Cabinets


Hanson cabinets are supplied in a variety of styles, sizes and options to meet your specific applications.

- Standing height base cabinets
- ADA height base cabinets
- Sitting height base cabinets
- Suspended cabinets
- Mobile cabinets
- Specialty cabinets
- Cabinet options and accessories



### Drawers

- Full access drawers feature a one piece drawer bodies for easy cleaning
- 18 ga. double walled on drawer body, with sound deadening full channeled edges
- Suspended on full extension, ball bearing, 100 lb. capacity drawer slides




# BASE CABINETS




Standing Height • 35-3/4"

hansonlab.com • 805-498-3121




## The HLS Advantages:

- Inset steel - flush face construction
- Standing height cabinets are 35-3/4" high and 22" deep
- Open cabinets have 2 adjustable shelves
- Cabinets with doors include one adjustable shelf
- Open and cabinets with doors have removable back panels
- See end of section for available accessories




Model	Width
AA15H	15"
AA18H	18-1/8"
AA24H	24-1/8"




Model	Width
AA15HL	15"
AA18HL	18-1/8"
AA24HL	24-1/8"

L= Left Hinged Door




Model	Width
AB15H	15"
AB18H	18-1/8"
AB24H	24-1/8"



Model	Width
AB15HL	15"
AB18HL	18-1/8"
AB24HL	24-1/8"

L= Left Hinged Door




Model	Width
AE15H	15"
AE18H	18-1/8"
AE24H	24-1/8"




Model	Width
AF18H	18-1/8"
AF24H	24-1/8"




Model	Width
AH18H	18-1/8"
AH24H	24-1/8"



Model	Width
AI18H	18-1/8"
AI24H	24-1/8"



Model	Width
AK18H	18-1/8"
AK24H	24-1/8"




Model	Width
AJ15H	15"
AJ18H	18-1/8"
AJ24H	24-1/8"




# BASE CABINETS

Standing Height • 35-3/4"




## The HLS Advantages:


- Inset steel - flush face construction
- Standing height cabinets are 35-3/4" high and 22" deep
- Open cabinets have 2 adjustable shelves
- Cabinets with doors include one adjustable shelf
- Open and cabinets with doors have removable back panels
- See end of section for available accessories




Model	Width
CE18H	18-1/8"
CE24H	24-1/8"




Model	Width
CF18H	18-1/8"
CF24H	24-1/8"




Model	Width
CA30H	30"
CA35H	35"
CA42H	42"
CA47H	47"



Model	Width
CD30H	30"
CD35H	35"
CD42H	42"
CD47H	47"




Model	Width
CG30H	30"
CG35H	35"
CG42H	42"
CG47H	47"




Model	Width
CG30HLC	30"
CG35HLC	35"
CG42HLC	42"
CG47HLC	47"


LC = Left Door



Model	Width
CN30H	30"
CN35H	35"
CN42H	42"
CN47H	47"



Model	Width
CP30H	30"
CP35H	35"
CP42H	42"
CP47H	47"



Model	Width
CT30H	30"
CT35H	35"
CT42H	42"
CT47H	47"




# BASE CABINETS




Standing Height • 35-3/4"

hansonlab.com • 805-498-3121




## The HLS Advantages:


- Inset steel - flush face construction
- Standing height cabinets are 35-3/4" high and 22" deep
- Open cabinets have 2 adjustable shelves
- Cabinets with doors include one adjustable shelf
- Open and cabinets with doors have removable back panels
- See end of section for available accessories




Model	Width
DY30H	30"
DY35H	35"
DY42H	42"
DY47H	47"




Model	Width
DA30H	30"
DA35H	35"
DA42H	42"
DA47H	47"



Model	Width
DD30H	30"
DD35H	35"
DD42H	42"
DD47H	47"




Model	Width
DC30H	30"
DC35H	35"
DC42H	42"
DC47H	47"




Model	Width
DE30H	30"
DE35H	35"
DE42H	42"
DE47H	47"



Model	Width
DF30H	30"
DF35H	35"
DF42H	42"
DF47H	47"




Model	Width
DG30H	30"
DG35H	35"
DG42H	42"
DG47H	47"




Model	Width
DG30HLC	30"
DG35HLC	35"
DG42HLC	42"
DG47HLC	47"

LC = Left Door




Model	Width
DH30H	30"
DH35H	35"
DH42H	42"
DH47H	47"




Model	Width
DH30HLC	30"
DH35HLC	35"
DH42HLC	42"
DH47HLC	47"

LC = Left Door



Model	Width
DI30H	30"
DI35H	35"
DI42H	42"
DI47H	47"




Model	Width
DI30HLC	30"
DI35HLC	35"
DI42HLC	42"
DI47HLC	47"

LC = Left Door


# BASE CABINETS

Standing Height • 35-3/4"




## The HLS Advantages:

- Inset steel - flush face construction
- Standing height cabinets are 35-3/4" high and 22" deep
- Open cabinets have 2 adjustable shelves
- Cabinets with doors include one adjustable shelf
- Open and cabinets with doors have removable back panels
- See end of section for available accessories




Model	Width
DK30H	30"
DK35H	35"
DK42H	42"
DK47H	47"




Model	Width
DK30HLC	30"
DK35HLC	35"
DK42HLC	42"
DK47HLC	47"

LC = Left Door




Model	Width
DO30H	30"
DO35H	35"
DO42H	42"
DO47H	47"




Model	Width
DO30HLC	30"
DO35HLC	35"
DO42HLC	42"
DO47HLC	47"


LC = Left Door



Model	Width
DN30H	30"
DN35H	35"
DN42H	42"
DN47H	47"



Model	Width
DP30H	30"
DP35H	35"
DP42H	42"
DP47H	47"




Model	Width
DQ30H	30"
DQ35H	35"
DQ42H	42"
DQ47H	47"




Model	Width
DR30H	30"
DR35H	35"
DR42H	42"
DR47H	47"



Model	Width
DS30H	30"
DS35H	35"
DS42H	42"
DS47H	47"



Model	Width
DT30H	30"
DT35H	35"
DT42H	42"
DT47H	47"



Model	Width
DZ30H	30"
DZ35H	35"
DZ42H	42"
DZ47H	47"



Model	Width
DZ30HLC	30"
DZ35HLC	35"
DZ42HLC	42"
DZ47HLC	47"



LC = Left Door

# CORNER BASE CABINETS




Standing Height • 35-3/4"

[hansonlab.com](http://hansonlab.com) • 805-498-3121




## The HLS Advantages:

- Inset steel - flush face construction
- Standing height cabinets are 35-3/4" high
- DX diagonal corner cabinets do not have shelves
- See end of section for available accessories




Model	Width
DW30H	30"
DW35H	35"
DW42H	42"
DW47H	47"



Model	Width
DW30HLC	30"
DW35HLC	35"
DW42HLC	42"
DW47HLC	47"


LC= Left Door



Model	Width
DX14H	32"


# SINK CABINETS

Standing Height • 35-3/4" - Sink




## The HLS Advantages:

- Inset steel - flush face construction
- Standing height cabinets are 35-3/4" high
- Sink cabinets have a lower rear panel
- If a cutout is required at the front toe kick for a foot-controlled fixture application, specify -FCCLP after the sink cabinet model number
- See end of section for available accessories




Model	Width
AC18H	18-1/8"
AC24H	24-1/8"




Model	Width
AC18HL	18-1/8"
AC24HL	24-1/8"

L = Left Hinged Door




Model	Width
DB30H	30"
DB35H	35"
DB42H	42"
DB47H	47"
DB52H	52"
DB58H	58"



Model	Width
DB47H3S	47"
DB52H3S	52"
DB58H3S	58"


Sitting Height • 29-3/4"

[hansonlab.com](http://hansonlab.com) • 805-498-3121




## The HLS Advantages:

- Inset steel - flush face construction
- Sitting height cabinets are 29-3/4" high and 22" deep
- Open cabinets have 2 adjustable shelves
- Cabinets with doors include one adjustable shelf
- Open and cabinets with doors have removable back panels
- See end of section for available accessories




Model	Width
AA15LL	15"
AA18LL	18-1/8"
AA24LL	24-1/8"


L = Left Hinged Door



Model	Width
AA15L	15"
AA18L	18-1/8"
AA24L	24-1/8"



Model	Width
AB15L	15"
AB18L	18-1/8"
AB24L	24-1/8"



Model	Width
AB15LL	15"
AB18LL	18-1/8"
AB24LL	24-1/8"

L = Left Hinged Door



Model	Width
AD15L	15"
AD18L	18-1/8"
AD24L	24-1/8"




Model	Width
AF18L	18-1/8"
AF24L	24-1/8"



Model	Width
AH18L	18-1/8"
AH24L	24-1/8"




Model	Width
AG18L	18-1/8"
AG24L	24-1/8"



Model	Width
AJ15L	15"
AJ18L	18-1/8"
AJ24L	24-1/8"


# BASE CABINETS

Sitting Height • 29-3/4"




## The HLS Advantages:


- Inset steel - flush face construction
- Sitting height cabinets are 29-3/4" high and 22" deep
- Open cabinets have 2 adjustable shelves
- Cabinets with doors include one adjustable shelf
- Open and cabinets with doors have removable back panels
- See end of section for available accessories




Model	Width
DY30L	30"
DY35L	35"
DY42L	42"
DY47L	47"




Model	Width
DA30L	30"
DA35L	35"
DA42L	42"
DA47L	47"




Model	Width
DD30L	30"
DD35L	35"
DD42L	42"
DD47L	47"




Model	Width
DC30L	30"
DC35L	35"
DC42L	42"
DC47L	47"




Model	Width
DE30L	30"
DE35L	35"
DE42L	42"
DE47L	47"



Model	Width
DF30L	30"
DF35L	35"
DF42L	42"
DF47L	47"




Model	Width
DG30L	30"
DG35L	35"
DG42L	42"
DG47L	47"




Model	Width
DG30LLC	30"
DG35LLC	35"
DG42LLC	42"
DG47LLC	47"

LC = Left Door




Model	Width
DH30L	30"
DH35L	35"
DH42L	42"
DH47L	47"




Model	Width
DH30LLC	30"
DH35LLC	35"
DH42LLC	42"
DH47LLC	47"

LC = Left Door



Model	Width
DI30L	30"
DI35L	35"
DI42L	42"
DI47L	47"



Model	Width
DI30LLC	30"
DI35LLC	35"
DI42LLC	42"
DI47LLC	47"


LC = Left Door

# BASE CABINETS




Sitting Height • 29-3/4"

hansonlab.com • 805-498-3121




## The HLS Advantages:

- Inset steel - flush face construction
- Sitting height cabinets are 29-3/4" high and 22" deep
- Open cabinets have 2 adjustable shelves
- Cabinets with doors include one adjustable shelf
- Open and cabinets with doors have removable back panels
- See end of section for available accessories




Model	Width
DK30L	30"
DK35L	35"
DK42L	42"
DK47L	47"




Model	Width
DK30LLC	30"
DK35LLC	35"
DK42LLC	42"
DK47LLC	47"

LC = Left Door




Model	Width
DO30L	30"
DO35L	35"
DO42L	42"
DO47L	47"




Model	Width
DO30LLC	30"
DO35LLC	35"
DO42LLC	42"
DO47LLC	47"


LC = Left Door




Model	Width
DN30L	30"
DN35L	35"
DN42L	42"
DN47L	47"




Model	Width
DP30L	30"
DP35L	35"
DP42L	42"
DP47L	47"




Model	Width
DQ30L	30"
DQ35L	35"
DQ42L	42"
DQ47L	47"



Model	Width
DR30L	30"
DR35L	35"
DR42L	42"
DR47L	47"




Model	Width
DJ30L	30"
DJ35L	35"
DJ42L	42"
DJ47L	47"




Model	Width
DJ30LLC	30"
DJ35LLC	35"
DJ42LLC	42"
DJ47LLC	47"

LC = Left Door



Model	Width
DZ30L	30"
DZ35L	35"
DZ42L	42"
DZ47L	47"




Model	Width
DZ30LLC	30"
DZ35LLC	35"
DZ42LLC	42"
DZ47LLC	47"

LC = Left Door


# BASE CABINETS

Sitting Height • 29-3/4"




## The HLS Advantages:


- Inset steel - flush face construction
- Sitting height cabinets are 29-3/4" high and 22" deep
- Open cabinets have 2 adjustable shelves
- Cabinets with doors include one adjustable shelf
- Open and cabinets with doors have removable back panels
- See end of section for available accessories



Model	Width
DS30L	30"
DS35L	35"
DS42L	42"
DS47L	47"



Model	Width
DT30L	30"
DT35L	35"
DT42L	42"
DT47L	47"




Model	Width
DV30L	30"
DV35L	35"
DV42L	42"
DV47L	47"

# CORNER BASE CABINETS



Sitting Height • 29-3/4"


[hansonlab.com](http://hansonlab.com) • 805-498-3121



## The HLS Advantages:


- Inset steel - flush face construction
- Sitting height cabinets are 29-3/4" high
- DX diagonal corner cabinets do not have shelves
- See end of section for available accessories

CORNER CABINET



Model	Width
DW30L	30"
DW35L	35"
DW42L	42"
DW47L	47"


CORNER CABINET



Model	Width
DW30LLC	30"
DW35LLC	35"
DW42LLC	42"
DW47LLC	47"

LC= Left Door



DIAGONAL CORNER CABINET



Model	Width
DX14L	32"


# SINK CABINETS

Sitting Height • 29-3/4"




## The HLS Advantages:

- Inset steel - flush face construction
- Sitting height cabinets are 29-3/4" high
- Sink cabinets have a lower rear panel
- If a cutout is required at the front toe kick for a foot-controlled fixture application, specify -FCCLP after the sink cabinet model number
- See end of section for available accessories




Model	Width
AC18L	18-1/8"
AC24L	24-1/8"




Model	Width
AC18LL	18-1/8"
AC24LL	24-1/8"


L = Left Hinged Door




Model	Width
DB30L	30"
DB35L	35"
DB42L	42"
DB47L	47"
DB52L	52"
DB58L	58"




Model	Width
DB47L3S	47"
DB52L3S	52"
DB58L3S	58"






### The HLS Advantages:


- Inset steel - flush face construction
- ADA height cabinets are 32" high and 22" deep
- Open cabinets have 2 adjustable shelves
- Cabinets with doors include one adjustable shelf
- Open and cabinets with doors have removable back panels
- See end of section for available accessories



Model	Width
AJ15ADA	15"
AJ18ADA	18-1/8"
AJ24ADA	24-1/8"




Model	Width
AA15ADA	15"
AA18ADA	18-1/8"
AA24ADA	24-1/8"




Model	Width
AA15ADAL	15"
AA18ADAL	18-1/8"
AA24ADAL	24-1/8"

L = Left Hinged Door



Model	Width
AB15ADA	15"
AB18ADA	18-1/8"
AB24ADA	24-1/8"



Model	Width
AB15ADAL	15"
AB18ADAL	18-1/8"
AB24ADAL	24-1/8"

L = Left Hinged Door



Model	Width
AF15ADA	15"
AF18ADA	18-1/8"
AF24ADA	24-1/8"




Model	Width
AG18ADA	18-1/8"
AG24ADA	24-1/8"




# BASE CABINETS

ADA Height • 32"




## The HLS Advantages:


- Inset steel - flush face construction
- ADA height cabinets are 32" high and 22" deep
- Open cabinets have 2 adjustable shelves
- Open and cabinets with doors have removable back panels
- See end of section for available accessories




Model	Width
DY30ADA	30"
DY35ADA	35"
DY42ADA	42"
DY47ADA	47"




Model	Width
DA30ADA	30"
DA35ADA	35"
DA42ADA	42"
DA47ADA	47"




Model	Width
DC30ADA	30"
DC35ADA	35"
DC42ADA	42"
DC47ADA	47"



Model	Width
DD30ADA	30"
DD35ADA	35"
DD42ADA	42"
DD47ADA	47"




Model	Width
DG30ADA	30"
DG35ADA	35"
DG42ADA	42"
DG47ADA	47"




Model	Width
DG30ADALC	30"
DG35ADALC	35"
DG42ADALC	42"
DG47ADALC	47"


LC = Left Door




Model	Width
DP30ADA	30"
DP35ADA	35"
DP42ADA	42"
DP47ADA	47"




Model	Width
DR30ADA	30"
DR35ADA	35"
DR42ADA	42"
DR47ADA	47"



Model	Width
DT30ADA	30"
DT35ADA	35"
DT42ADA	42"
DT47ADA	47"



Model	Width
DU30ADA	30"
DU35ADA	35"
DU42ADA	42"
DU47ADA	47"



Model	Width
DU30ADALC	30"
DU35ADALC	35"
DU42ADALC	42"
DU47ADALC	47"


LC = Left Door

# CORNER BASE CABINETS




ADA Height • 32"

[hansonlab.com](http://hansonlab.com) • 805-498-3121




## The HLS Advantages:

- Inset steel - flush face construction
- ADA height cabinets are 32" high
- Corner cabinets have one adjustable shelf




Model	Width
DW30ADA	30"
DW35ADA	35"
DW42ADA	42"
DW47ADA	47"




Model	Width
DW30ADALC	30"
DW35ADALC	35"
DW42ADALC	42"
DW47ADALC	47"

LC = Left Door



# SINK CABINETS


ADA Height • 32"



## The HLS Advantages:


- Inset steel - flush face construction
- ADA height cabinets are 32" high
- Sink cabinets without a bottom floor have a vented integral toe-kick for wheelchair accessibility
- Sink cabinet style with floor are not wheelchair accessible
- Sink cabinets have a lower rear panel, no shelves

WITHOUT FLOOR  
WHEELCHAIR ACCESSIBLE



Model	Width
DB35ADA	35"
DB42ADA	42"
DB47ADA	47"
DB52ADA	52"
DB58ADA	58"


WITH FLOOR  
NON WHEELCHAIR ACCESSIBLE



Model	Width
AC18ADAF	18-1/8"
AC24ADAF	24-1/8"

F = Floor


WITH FLOOR  
NON WHEELCHAIR ACCESSIBLE



Model	Width
AC18ADALF	18-1/8"
AC24ADALF	24-1/8"

L = Left Hinged Door  
F = Floor

WITH FLOOR  
NON WHEELCHAIR ACCESSIBLE





Model	Width
DB30ADAF	30"
DB35ADAF	35"
DB42ADAF	42"
DB47ADAF	47"
DB52ADAF	52"
DB58ADAF	58"

F = Floor


## The HLS Advantages

- Inset steel - flush face construction
- For use with HLS Agility Adaptable, M-Series workstations & HDW tables\*
- Suspended cabinets are 21" high and 22" deep
- Open cabinets and cabinets with doors have one adjustable shelf
- Cabinets can be removed, relocated, or moved to any location on the table chosen
- Suspended cabinets have finished backs and include cabinet attachment brackets

\*specify suspended cabinet option when ordering workstation / table



Model	Width
AJ15U	15"
AJ18U	18-1/8"
AJ24U	24-1/8"




Model	Width
AA15U	15"
AA18U	18-1/8"
AA24U	24-1/8"




Model	Width
AA15UL	15"
AA18UL	18-1/8"
AA24UL	24-1/8"


L = Left Hinged Door



Model	Width
DA30U	30"
DA35U	35"
DA42U	42"




Model	Width
AB15U	15"
AB18U	18-1/8"
AB24U	24-1/8"




Model	Width
AB15UL	15"
AB18UL	18-1/8"
AB24UL	24-1/8"


L = Left Hinged Door



Model	Width
AF18U	18-1/8"
AF24U	24-1/8"



Model	Width
AG18U	18-1/8"
AG24U	24-1/8"




Model	Width
DD30U	30"
DD35U	35"
DD42U	42"

# MOBILE CABINETS


Standard Height • 32 1/2"

## The HLS Advantages

- Inset steel - flush face construction
- Painted steel top and finished backs are standard
- Mobile height cabinets are 32-1/2" high and 22" deep
- Open cabinets have 2 adjustable shelves
- Cabinets with doors have one adjustable shelf
- 3" dia. wheels on 4-1/4" high casters (front casters swivel w/ brake; rear casters fixed, no brake)
- Optional counter tops available - specify material type and thickness (add thickness to mobile cabinet height)




Model	Width
AA15HM	15"
AA18HM	18-1/8"
AA24HM	24-1/8"




Model	Width
AA15HML	15"
AA18HML	18-1/8"
AA24HML	24-1/8"

L = Left Hinged Door



Model	Width
AB15HM	15"
AB18HM	18-1/8"
AB24HM	24-1/8"




Model	Width
AB15HML	15"
AB18HML	18-1/8"
AB24HML	24-1/8"


L = Left Hinged Door




Model	Width
AF15HM	15"
AF18HM	18-1/8"
AF24HM	24-1/8"




Model	Width
AG18HM	18"
AG24HM	24-1/8"




Model	Width
AJ15HM	15"
AJ18HM	18-1/8"
AJ24HM	24-1/8"



Model	Width
CG30HM	30"
CG35HM	35"
CG42HM	42"
CG47HM	47"



Model	Width
CG30HMLC	30"
CG35HMLC	35"
CG42HMLC	42"
CG47HMLC	47"




Model	Width
DA30HM	30"
DA35HM	35"
DA42HM	42"
DA47HM	47"


LC = Left Door




Model	Width
DD30HM	30"
DD35HM	35"
DD42HM	42"
DD47HM	47"



Model	Width
DG30HM	30"
DG35HM	35"
DG42HM	42"
DG47HM	47"



Model	Width
DG30HMLC	30"
DG35HMLC	35"
DG42HMLC	42"
DG47HMLC	47"





Model	Width
DY30HM	30"
DY35HM	35"
DY42HM	42"
DY47HM	47"


LC = Left Door

## The HLS Advantages


- Inset steel - flush face construction
- Painted steel top and finished backs are standard
- Mobile height cabinets are 25-1/4" high and 22" deep
- Open cabinets have 2 adjustable shelves
- Cabinets with doors have one adjustable shelf
- 3" dia. wheels on 4-1/4" high casters (front casters swivel w/ brake; rear casters fixed, no brake)
- Optional counter tops available - specify material type and thickness (add thickness to mobile cabinet height)



Model	Width
AJ15LM	15"
AJ18LM	18-1/8"
AJ24LM	24-1/8"




Model	Width
AA15LM	15"
AA18LM	18-1/8"
AA24LM	24-1/8"



Model	Width
AA15LML	15"
AA18LML	18-1/8"
AA24LML	24-1/8"

L = Left Hinged Door



Model	Width
DA30LM	30"
DA35LM	35"
DA42LM	42"



Model	Width
AB15LM	15"
AB18LM	18-1/8"
AB24LM	24-1/8"




Model	Width
AB15LML	15"
AB18LML	18-1/8"
AB24LML	24-1/8"


L = Left Hinged Door



Model	Width
AF18LM	18-1/8"
AF24LM	24-1/8"



Model	Width
AG18LM	18-1/8"
AG24LM	24-1/8"




Model	Width
DD30LM	30"
DD35LM	35"
DD42LM	42"


# CHEMICAL STORAGE CABINETS

Standing Height • 35-3/4"


## Acid / Base / Corrosive Storage Cabinets



- Designed for under counter and under fume hood applications, for storage of corrosive materials
- Pneumatic self-closing door - doors have a 165°F fusible link that automatically closes the doors in case of a fire
- 3 point locking door mechanism with keyed lock
- Includes one adjustable shelf
- Two vent holes are provided in the back of cabinet for use with an optional VKACID vent kit
- Lined with corrosive resistant liner board
- Each exterior door face includes 9" x 2-1/4" adhesive label 'ACID' - located as shown
- Cabinet also available with 'BASE' label. Use BSB part number prefix in place of ASB




Model	Width
ASB18H	18"
ASB23H	23"
ASB24H	24"




Model	Width
ASB18HL	18"
ASB23HL	23"
ASB24HL	24"

L = Left Hinged Door



Model	Width
ASB30H	30"
ASB35H	35"
ASB47H	47"

Model	Capacity
ASB18H	10 Gal.
ASB23H	12 Gal.
ASB24H	12 Gal.
ASB30H	16 Gal.
ASB35H	18 Gal.
ASB47H	27 Gal.



## Acid / Corrosive Cabinet Vent Kit


- For connection from acid cabinet to fume hood
- Includes an 12' long, 1" diameter O.D. PVC flexible hose
- A connection is provided for attachment to the back of the cabinet

Note - these kits may be included with cabinet assemblies listed above

Model	Hose Length
VKACID	12'
VKACID-VAV	25'

Standing Height • 35-3/4"

hansonlab.com • 805-498-3121




## Flammable Storage Cabinets

- Use under counters, under fume hoods or in stand-alone applications
- Double-wall 18 ga. construction w / 1.5" air space
- Pneumatic self-closing door - doors have a 165°F fusible link that automatically closes the door in case of a fire
- 3 point locking door mechanism with keyed lock
- Chemical-resistant interior and exterior finish
- Includes one adjustable shelf
- 18 ga., 2" sill at bottom to retain spills
- Two spark arrester vent holes are provided in the back of cabinet for use with optional VKFLAM vent kits
- A removable panel is provided in the back of the cabinet for access to plumbing
- Each exterior door face includes 9" x 2-1/4" adhesive label 'FLAMMABLE' - located as shown




Model	Width
FSB18H	18"
FSB24H	24"




Model	Width
FSB18HL	18"
FSB24HL	24"

L = Left Hinged Door



Model	Width
FSB30H	30"
FSB35H	35"
FSB47H	47"

Model	Capacity
FSB18H	10 Gal.
FSB24H	12 Gal.
FSB30H	16 Gal.
FSB35H	18 Gal.
FSB47H	27 Gal.




## Flammable Cabinet Vent Kit (optional)

- For vent connection from flammable cabinet
- Includes a 12' long, 1" diameter O.D. galvanized steel flexible hose
- Connection by others
- A connection is provided for attachment to the back of the cabinet

Model	Hose Length
VKFLAM	12'
VKFLAM-VAV	25'


# CHEMICAL STORAGE CABINETS

Sitting Height • 29-3/4"




## Acid / Base / Corrosive Storage Cabinets

- Designed for under counter and under fume hood applications, for storage of corrosive materials
- Pneumatic self-closing door - doors have a 165°F fusible link that automatically closes the doors in case of a fire
- 3 point locking door mechanism with keyed lock
- Includes one adjustable shelf
- Two vent holes are provided in the back of cabinet for use with optional VKACID vent kit
- Lined with corrosive resistant liner board
- Each exterior door face includes 9" x 2-1/4" adhesive label 'ACID' - located as shown
- Cabinet also available with 'BASE' label. Use BSB part number prefix in place of ASB




Model	Width
ASB18L	18"
ASB23L	23"
ASB24L	24"




Model	Width
ASB18LL	18"
ASB23LL	23"
ASB24LL	24"

L = Left Hinged Door



Model	Width
ASB30L	30"
ASB35L	35"
ASB47L	47"

Model	Capacity
ASB18L	8 Gal.
ASB23L	10 Gal.
ASB24L	10 Gal.
ASB30L	12 Gal.
ASB35L	16 Gal.
ASB47L	25 Gal.




## Acid / Corrosive Cabinet Vent Kit

- For connection from acid cabinet to fume hood
- Includes a 12' long, 1" diameter O.D. PVC flexible hose
- A connection is provided for attachment to the back of cabinet

Model	Hose Length
VKACID	12'
VKACIDVAV	25'


Sitting Height • 29-3/4"

hansonlab.com • 805-498-3121




## Flammable Storage Cabinets

- Use under counters, under fume hoods or in stand-alone applications
- Double-wall 18 ga. construction w / 1.5" air space
- Pneumatic self-closing door - doors have a 165°F fusible link that automatically closes the door in case of a fire
- 3 point locking door mechanism with keyed lock
- Chemical-resistant interior and exterior finish
- Includes one adjustable shelf
- 18 ga., 2" sill at bottom to retain spills
- Two spark arrester vent holes are provided in the back of cabinet for use with optional VKFLAM vent kits
- A removable panel is provided in the back of the cabinet for access to plumbing
- Each exterior door face includes 9" x 2-1/4" adhesive label 'FLAMMABLE' - located as shown




Model	Width
FSB18L	18"
FSB24L	24"




Model	Width
FSB18LL	18"
FSB24LL	24"

L = Left Hinged Door



Model	Width
FSB30L	30"
FSB35L	35"
FSB47L	47"

Model	Capacity
FSB18L	8 Gal.
FSB24L	10 Gal.
FSB30L	12 Gal.
FSB35L	16 Gal.
FSB47L	25 Gal.




## Flammable Cabinet Vent Kit (optional)

- For vent connection from flammable cabinet
- Includes a 12' long, 1" diameter I.D. galvanized steel flexible hose
- Connection by others
- A connection is provided for attachment to the back of cabinet

Model	Hose Length
VKFLAM	12'
VKFLAM-VAV	25'


# CHEMICAL STORAGE CABINETS

ADA Height • 32"




## Acid / Base / Corrosive Storage Cabinets

- Designed for under counter and under fume hood applications, for storage of corrosive materials
- Pneumatic self-closing door - doors have a 165°F fusible link that automatically closes the doors in case of a fire
- 3 point locking door mechanism with keyed lock
- Includes one adjustable shelf
- Two vent holes are provided in the back of cabinet for use with optional VKACID vent kit
- Lined with corrosive resistant liner board
- Each exterior door face includes 9" x 2-1/4" adhesive label 'ACID' - located as shown
- Cabinet also available with 'BASE' label. Use BSB part number prefix in place of ASB




Model	Width
ASB18ADA	18"
ASB23ADA	23"
ASB24ADA	24"




Model	Width
ASB18ADAL	18"
ASB23ADAL	23"
ASB24ADAL	24"

L = Left Hinged Door



Model	Width
ASB30ADA	30"
ASB35ADA	35"
ASB47ADA	47"


Model	Capacity
ASB18ADA	8 Gal.
ASB23ADA	10 Gal.
ASB24ADA	10 Gal.
ASB30ADA	12 Gal.
ASB35ADA	16 Gal.
ASB47ADA	25 Gal.



## Acid / Corrosive Cabinet Vent Kit


- For connection from acid cabinet to fume hood
- Includes a 12' long, 1" diameter O.D. PVC flexible hose
- A connection is provided for attachment to the back of cabinet

Model	Hose Length
VKACID	12'
VKACID-VAV	25'




## Flammable Storage Cabinets

- Use under counters, under fume hoods or in stand-alone applications
- Double-wall 18 ga. construction w / 1.5" air space
- Pneumatic self-closing door - doors have a 165°F fusible link that automatically closes the door in case of a fire
- 3 point locking door mechanism with keyed lock
- Chemical-resistant interior and exterior finish
- Includes one adjustable shelf
- 18 ga., 2" sill at bottom to retain spills
- Two spark arrester vent holes are provided in the back of cabinet for use with optional VKFLAM vent kits
- A removable panel is provided in the back of the cabinet for access to plumbing
- Each exterior door face includes 9" x 2-1/4" adhesive label 'FLAMMABLE' - located as shown




Model	Width
FSB18ADA	18"
FSB24ADA	24"




Model	Width
FSB18ADAL	18"
FSB24ADAL	24"

L = Left Hinged Door



Model	Width
FSB30ADA	30"
FSB35ADA	35"
FSB47ADA	47"


Model	Capacity
FSB18ADA	8 Gal.
FSB24ADA	10 Gal.
FSB30ADA	12 Gal.
FSB35ADA	16 Gal.
FSB47ADA	25 Gal.



## Acid / Corrosive Cabinet Vent Kit

- For connection from acid cabinet to fume hood
- Includes a 12' long, 1" diameter O.D. PVC flexible hose
- A connection is provided for attachment to the back of cabinet

Model	Hose Length
VKACID	12'
VKACID-VAV	25'




## Flammable Cabinet Vent Kit (optional)

- For vent connection from flammable cabinet
- Includes a 12' long, 1" diameter I.D. galvanized steel flexible hose
- Connection by others
- A connection is provided for attachment to the back of cabinet

Model	Hose Length
VKFLAM	12'
VKFLAM-VAV	25'


# SPECIALTY CABINETS

Standing Height • 35-3/4"




## Equipment Cabinet

- Used for large or roll-in equipment
- Cabinet has a integral vented toe kick without a floor
- Full height back with removable back panel
- Manufactured in 22" and 28" deep cabinets




Model	Depth	Width
CU47HA	22"	47"
CU47HB	28"	47"
CU52HA	22"	52"
CU52HB	28"	52"
CU58HA	22"	58"
CU58HB	28"	58"




## Vacuum Pump Cabinets

- Models with the suffix "F" come with a floor and pull out bottom shelf
- Designed with pre-wired electrical package, illuminated light switch, interior duplex receptacle, fan, and vented toekick (35", 42", and 47" wide cabinets have 1 electrical package provided with 2 fans)
- One 15ft. electrical conduit whip is provided with each cabinet
- Sound deadening interior insulation
- A vent hole and removable access panel are provided in the back of cabinet for use with included VKVAC vent kit
- See page 42 for optional mobile platform




Model	Width
VAC23HF	23"
VAC24HF	24"




Model	Width
VAC23HLF	23"
VAC24HLF	24"

L = Left Hinged Door





Model	Width
VAC30HF	30"
VAC35HF	35"
VAC42HF	42"
VAC47HF	47"




Model	Width
VAC30HF-R	30"
VAC35HF-R	35"
VAC42HF-R	42"
VAC47HF-R	47"

R = Right Hand Switch




Model	Width
VAC23H	23"
VAC24H	24"




Model	Width
VAC23HL	23"
VAC24HL	24"

L = Left Hinged Door



Model	Width
VAC30H	30"
VAC35H	35"
VAC42H	42"
VAC47H	47"




Model	Width
VAC30H-R	30"
VAC35H-R	35"
VAC42H-R	42"
VAC47H-R	47"

R = Right Hand Switch


# SPECIALTY CABINETS

Sitting Height • 29-3/4"




## Vacuum Pump Cabinets

- Models with the suffix "F" come with a floor and pull out bottom shelf
- Designed with pre-wired electrical package, illuminated light switch, interior duplex receptacle, fan, and vented toekick (35", 42", and 47" wide cabinets have 1 electrical package provided with 2 fans)
- One 15ft. electrical conduit whip is provided with each cabinet
- Sound deadening interior insulation
- A vent hole and removable access panel are provided in the back of cabinet for use with included VKVAC vent kit
- See page 42 for optional mobile platform




Model	Width
VAC23LF	23"
VAC24LF	24"




Model	Width
VAC23LLF	23"
VAC24LLF	24"

L = Left Hinged Door





Model	Width
VAC30LF	30"
VAC35LF	35"
VAC42LF	42"
VAC47LF	47"




Model	Width
VAC30LF-R	30"
VAC35LF-R	35"
VAC42LF-R	42"
VAC47LF-R	47"

R = Right Hand Switch




Model	Width
VAC23L	23"
VAC24L	24"




Model	Width
VAC23LL	23"
VAC24LL	24"

L = Left Hinged Door




Model	Width
VAC30L	30"
VAC35L	35"
VAC42L	42"
VAC47L	47"




Model	Width
VAC30L-R	30"
VAC35L-R	35"
VAC42L-R	42"
VAC47L-R	47"

R = Right Hand Switch




## Vacuum Pump Cabinets

- Models with the suffix "F" come with a floor and pull out bottom shelf
- Designed with pre-wired electrical package, illuminated light switch, interior duplex receptacle, fan, and vented toekick (35", 42", and 47" wide cabinets have 1 electrical package provided with 2 fans)
- One 15 ft. electrical conduit whip is provided with each cabinet
- Sound deadening interior insulation
- A vent hole and removable access panel are provided in the back of cabinet for use with included VKVAC vent kit
- See page 42 for optional mobile platform




Model	Width
VAC23ADAF	23"
VAC24ADAF	24"




Model	Width
VAC23ADALF	18"
VAC24ADALF	23"

L = Left Hinged Door





Model	Width
VAC30ADAF	30"
VAC35ADAF	35"
VAC42ADAF	42"
VAC47ADAF	47"




Model	Width
VAC30ADAF-R	30"
VAC35ADAF-R	35"
VAC42ADAF-R	42"
VAC47ADAF-R	47"

R = Right Hand Switch




Model	Width
VAC23ADA	23"
VAC24ADA	24"




Model	Width
VAC23ADAL	23"
VAC24ADAL	24"

L = Left Hinged Door



Model	Width
VAC30ADA	30"
VAC35ADA	35"
VAC42ADA	42"
VAC47ADA	47"



Model	Width
VAC30ADA-R	30"
VAC35ADA-R	35"
VAC42ADA-R	42"
VAC47ADA-R	47"


R = Right Hand Switch

# CASEWORK ACCESSORIES

## Cabinet Options

### Cabinet Locks


- Used on drawers or doors for base, wall, and tall cabinets (not available on ADA sink cabinets)
- Specify -LK at the end of the cabinet model number for locks on all doors and drawers on the cabinet (please specify at time of order if locks are to be **keyed alike**, or **keyed different**)
- Not available on ADA sink cabinets



Model  
-LK

### Finished Backs


- Used to finish the back of base cabinets
- Installed at the factory
- Specify -FB at the end of the cabinet model number



Model  
-FB

### Pull Out Shelves


- Install in base cabinet with open doors or open cabinets
- Pull out shelves have 2" high lips on all 4 sides
- Unless otherwise specified, the shelf will be at bottom of cabinet
- Specify -PS at the end of the cabinet model number



Model  
-PS

### Pull Boards


- Pull boards are available on most base cabinets
- Not available on cabinets that do not have a full height center vertical rail
- Specify for left, right, or both sides for 30" through 47' wide base cabinets
- To add pull boards, specify -PB at the end of the model number



Model  
-PB  
-PBL(Left)  
-PBR(Right)  
-PBB(Both)

### Mobile Platform


- Mobile platform for use with vacuum pump cabinets
- Platform measures 16" (front-to-rear) with 5-1/2" overall height and 2" tall vertical retaining lips on 4 sides
- Includes four casters



Model	Width
MP-23	18" w Fits 23" Cabinet
MP-24	19" w Fits 24" Cabinet
MP-30	25" w Fits 30" Cabinet
MP-35	30" w Fits 35" Cabinet
MP-42	37" w Fits 42" Cabinet
MP-47	42" w Fits 47" Cabinet

### Lazy Susan

- Fits into DX style corner cabinet




Model	Description
LS-1	1-Tier "Lazy Susan"
LS-2	2 Tier "Lazy Susan"

### Base Filler Panel


- Used between 2 cabinets or between a cabinet and the wall
  - Specify -T at the end of the number if the panel needs to be field trimmed with an SA1 double-slip and single angle
- NOTE: 1" wide fillers are not available with -T variation

#### Standing Height Filler Panel




Model	Width	Height
BF1M	1"	35-3/4"
BF4M	4"	35-3/4"
BF6M	6"	35-3/4"
BF8M	8"	35-3/4"
BF12M	12"	35-3/4"
BF14M	14"	35-3/4"
BF18M	18"	35-3/4"

#### ADA Height Filler Panel




Model	Width	Height
BF1ADA	1"	32"
BF4ADA	4"	32"
BF6ADA	6"	32"
BF8ADA	8"	32"
BF12ADA	12"	32"
BF14ADA	14"	32"
BF18ADA	18"	32"

#### Low Height Filler Panel




Model	Width	Height
BF1L	1"	29-3/4"
BF4L	4"	29-3/4"
BF6L	6"	29-3/4"
BF8L	8"	29-3/4"
BF12L	12"	29-3/4"
BF14L	14"	29-3/4"
BF18L	18"	29-3/4"

### End Filler Panel



Model	Width	Height
PSF1L	1"	29-3/4"
PSF7L	7"	29-3/4"
PSF13L	13"	29-3/4"
PSF1ADA	1"	32"
PSF7ADA	7"	32"
PSF13ADA	13"	32"
PSF1M	1"	35-3/4"
PSF7M	7"	35-3/4"
PSF13M	13"	35-3/4"


### Island Filler Panel



Model	Width	Height
BCF2L	2"	29-3/4"
BCF7L	7"	29-3/4"
BCF14L	14"	29-3/4"
BCF26L	26"	29-3/4"
BCF2ADA	2"	32"
BCF7ADA	7"	32"
BCF14ADA	14"	32"
BCF26ADA	26"	32"
BCF2M	2"	35-3/4"
BCF7M	7"	35-3/4"
BCF14M	14"	35-3/4"
BCF26M	26"	35-3/4"

### Double-Slip Angle

- Field trimmable for use with any height cabinet



Model	Width	Height
SA1L	1"	25-3/4"
SA1ADA	1"	28-1/8"
SA1	1"	31-3/4"


Includes a single angle.

# CASEWORK ACCESSORIES

## Knee Space Panels

### KSP Filler Panel - applications w/o Apron Rail

- Full height panel, used when there is no rear apron or drawer frame
- For use between 2 cabinets or between a cabinet and wall
- Includes two 1" x 1" angles for installation
- KSPs can be field trimmed




Model	Width	Height
KSP1829	18"	29"
KSP2429	24"	29"
KSP3029	30"	29"
KSP3529	35"	29"
KSP4229	42"	29"
KSP4729	47"	29"
KSP5829	58"	29"
KSP7229	72"	29"

Model	Width	Height
KSP1832	18"	32"
KSP2432	24"	32"
KSP3032	30"	32"
KSP3532	35"	32"
KSP4232	42"	32"
KSP4732	47"	32"
KSP5832	58"	32"
KSP7232	72"	32"

Model	Width	Height
KSP1835	18"	35"
KSP2435	24"	35"
KSP3035	30"	35"
KSP3535	35"	35"
KSP4235	42"	35"
KSP4735	47"	35"
KSP5835	58"	35"
KSP7235	72"	35"

### KSP Filler Panel - applications w/ Apron Rail

- Used when a rear apron or drawer frame is used
- For use between 2 cabinets or between a cabinet and wall
- Includes two 1" x 1" angles for installation
- KSPs can be field trimmed



#### Use w/ 29"H Cabinets

Model	Width	Height
KSP1826	18"	26"
KSP2426	24"	26"
KSP3026	30"	26"
KSP3526	35"	26"
KSP4226	42"	26"
KSP4726	47"	26"
KSP5826	58"	26"
KSP7226	72"	26"

#### Use w/ 32"H Cabinets


Model	Width	Height
KSP1828	18"	28"
KSP2428	24"	28"
KSP3028	30"	28"
KSP3528	35"	28"
KSP4228	42"	28"
KSP4728	47"	28"
KSP5828	58"	28"
KSP7228	72"	28"

#### Use w/ 35"H Cabinets

Model	Width	Height
KSP1832	18"	32"
KSP2432	24"	32"
KSP3032	30"	32"
KSP3532	35"	32"
KSP4232	42"	32"
KSP4732	47"	32"
KSP5832	58"	32"
KSP7232	72"	32"

### KSP Filler Panel - applications w/o Apron Rail

- Full height panel, used when there is no rear apron or drawer frame
- For use between 2 cabinets or between a cabinet and wall
- Includes two 1" x 1" angles for installation
- KSPs can be field trimmed




Model	Width	Height
KSP1829E	18"	29"
KSP2429E	24"	29"
KSP3029E	30"	29"
KSP3529E	35"	29"
KSP4229E	42"	29"
KSP4729E	47"	29"
KSP5829E	58"	29"
KSP7229E	72"	29"

Model	Width	Height
KSP1832E	18"	32"
KSP2432E	24"	32"
KSP3032E	30"	32"
KSP3532E	35"	32"
KSP4232E	42"	32"
KSP4732E	47"	32"
KSP5832E	58"	32"
KSP7232E	72"	32"

Model	Width	Height
KSP1835E	18"	35"
KSP2435E	24"	35"
KSP3035E	30"	35"
KSP3535E	35"	35"
KSP4235E	42"	35"
KSP4735E	47"	35"
KSP5835E	58"	35"
KSP7235E	72"	35"

### KSP Filler Panel - applications w/ Apron Rail

- Used when a rear apron or drawer frame is used
- For use between 2 cabinets or between a cabinet and wall
- Includes two 1" x 1" angles for installation
- KSPs can be field trimmed



#### Use w/ 29"H Cabinets

Model	Width	Height
KSP1826E	18"	26"
KSP2426E	24"	26"
KSP3026E	30"	26"
KSP3526E	35"	26"
KSP4226E	42"	26"
KSP4726E	47"	26"
KSP5826E	58"	26"
KSP7226E	72"	26"

#### Use w/ 32"H Cabinets

Model	Width	Height
KSP1828E	18"	28"
KSP2428E	24"	28"
KSP3028E	30"	28"
KSP3528E	35"	28"
KSP4228E	42"	28"
KSP4728E	47"	28"
KSP5828E	58"	28"
KSP7228E	72"	28"

#### Use w/ 35"H Cabinets


Model	Width	Height
KSP1832E	18"	32"
KSP2432E	24"	32"
KSP3032E	30"	32"
KSP3532E	35"	32"
KSP4232E	42"	32"
KSP4732E	47"	32"
KSP5832E	58"	32"
KSP7232E	72"	32"

# CASEWORK ACCESSORIES

## Rear Enclosure Panels – REP

### Rear Enclosure Panels

- Rear enclosure panels are used to enclose the back of cabinets or island benches when the counter depth is not deep enough for back-to-back cabinets
- Panels are manufactured in Standing, ADA, and Sitting height panels to match our cabinet heights and include leveling bolts in each corner
- For 2 sided toekick specify 2S at the end of the product number



### Sitting Height 29-3/4" x 6" Deep

Model	Width	Depth	Height
REP18L6	18-1/8"	6"	29-3/4"
REP24L6	24-1/8"	6"	29-3/4"
REP30L6	30"	6"	29-3/4"
REP35L6	35"	6"	29-3/4"
REP42L6	42"	6"	29-3/4"
REP47L6	47"	6"	29-3/4"
REP58L6	58"	6"	29-3/4"

### Sitting Height 29-3/4" x 12" Deep

Model	Width	Depth	Height
REP18L12	18-1/8"	12"	29-3/4"
REP24L12	24-1/8"	12"	29-3/4"
REP30L12	30"	12"	29-3/4"
REP35L12	35"	12"	29-3/4"
REP42L12	42"	12"	29-3/4"
REP47L12	47"	12"	29-3/4"
REP58L12	58"	12"	29-3/4"

### ADA Height 32" x 6" Deep

Model	Width	Depth	Height
REP18A6	18-1/8"	6"	32"
REP24A6	24-1/8"	6"	32"
REP30A6	30"	6"	32"
REP35A6	35"	6"	32"
REP42A6	42"	6"	32"
REP47A6	47"	6"	32"
REP58A6	58"	6"	32"

### ADA Height 32" x 12" Deep

Model	Width	Depth	Height
REP18A12	18-1/8"	12"	32"
REP24A12	24-1/8"	12"	32"
REP30A12	30"	12"	32"
REP35A12	35"	12"	32"
REP42A12	42"	12"	32"
REP47A12	47"	12"	32"
REP58A12	58"	12"	32"


### Standing Height 35-3/4" x 6" Deep

Model	Width	Depth	Height
REP18H6	18-1/8"	6"	35-3/4"
REP24H6	24-1/8"	6"	35-3/4"
REP30H6	30"	6"	35-3/4"
REP35H6	35"	6"	35-3/4"
REP42H6	42"	6"	35-3/4"
REP47H6	47"	6"	35-3/4"
REP58H6	58"	6"	35-3/4"

### Standing Height 35-3/4" x 12" Deep


Model	Width	Depth	Height
REP18H12	18-1/8"	12"	35-3/4"
REP24H12	24-1/8"	12"	35-3/4"
REP30H12	30"	12"	35-3/4"
REP35H12	35"	12"	35-3/4"
REP42H12	42"	12"	35-3/4"
REP47H12	47"	12"	35-3/4"
REP58H12	58"	12"	35-3/4"

- 6" deep panels are typically used with 30" deep counter tops for a 1" front and rear overhang




TOP VIEW

- 12" deep panels are typically used with 36" deep counter tops for a 1" front and rear overhang






TOP VIEW



### Apron Rails and Drawer Frames

- Apron rails and drawer frames can be used as a knee space between two cabinets, between a cabinet and a wall, or between a cabinet and legs or end support panels
- Apron rails or drawer frames are NOT designed to be used between two sets of legs as a freestanding table. For freestanding tables see our Agility / Mobility section of the catalog
- Drawer frames have a 3-1/2" high frame with a 2" high drawer



#### Apron Rails

Model	Height	Width
A24	3-1/2"	24"
A30	3-1/2"	30"
A35	3-1/2"	35"
A42	3-1/2"	42"
A47	3-1/2"	47"
A54	3-1/2"	54"
A60	3-1/2"	60"
A66	3-1/2"	66"
A72	3-1/2"	72"
A224	2"	24"
A230	2"	30"
A235	2"	35"
A242	2"	42"
A247	2"	47"
A254	2"	54"
A260	2"	60"
A266	2"	66"
A272	2"	72"

#### Drawer Frames One Pencil Drawer

Model	Width	Depth
DH-1-24A	24"	22"
DH-1-24B	24"	28"
DH-1-30A	30"	22"
DH-1-30B	30"	28"
DH-1-35A	35"	22"
DH-1-35B	35"	28"
DH-1-42A	42"	22"
DH-1-42B	42"	28"
DH-1-47A	47"	22"
DH-1-47B	47"	28"


#### Drawer Frames Two Pencil Drawer

Model	Width	Depth
DH-2-47A	47"	22"
DH-2-47B	47"	28"
DH-2-54A	54"	22"
DH-2-54B	54"	28"
DH-2-58A	58"	22"
DH-2-58B	58"	28"
DH-2-60A	60"	22"
DH-2-60B	60"	28"
DH-2-66A	66"	22"
DH-2-66B	66"	28"

### Toe Space Filler Panel

- Used at corner where cabinets meet at a 90° angle
- 18 ga. steel painted to match cabinets


Model  
TSF



# CASEWORK ACCESSORIES

## Leg Assemblies & Countertop Supports

### Single Leg



- Used with apron rails or drawer frames
- 2" x 2" 16 ga. construction
- Leg levelers and vinyl shoe included

Model	Height
L30	29-3/4"
L32	32"
L36	35-3/4"


### Single Flanged Leg



- Used to provide added support of counter tops
- 4" x 4" plate, offset at the top with holes for attachment to counter top
- Leg levelers and vinyl shoe included

Model	Height
FL30	29-3/4"
FL32	32"
FL36	35-3/4"


### H-Leg



- Used with apron rails or drawer frames
- Leg levelers and vinyl shoes provided

Model	Height	Depth
HL30A	29-3/4"	22"
HL30B	29-3/4"	28"
HL32A	32"	22"
HL32B	32"	28"
HL36A	35-3/4"	22"
HL36B	35-3/4"	28"


### Flanged H-Leg



- Used to provide added support of counter tops
- 4" x 4" plates, offset at the top with holes for attachment to countertop
- Leg levelers and vinyl shoe included

Model	Height	Depth
FHL30A	29-3/4"	22"
FHL30B	29-3/4"	28"
FHL32A	32"	22"
FHL32B	32"	28"
FHL36A	35-3/4"	22"
FHL36B	35-3/4"	28"


### Countertop Support Brackets



- Painted to match casework
- In-wall backing required for support

Model	Height	Depth
CSB21	20"	21"
CSB27	20"	27"
EH1824	18"	24"

### End Support Panel



- Used in place of legs on one end of apron rails or drawer frames
- 1-1/4" wide
- Leg levelers provided


Model	Height	Depth
HGL30A	29-3/4"	22"
HGL30B	29-3/4"	28"
HGL32A	32"	22"
HGL32B	32"	28"
HGL36A	35-3/4"	22"
HGL36B	35-3/4"	28"

### Function & Design

- HLS offers a variety of Umbilical Pipe Enclosures used to conceal service lines that come down from the ceiling at either walls or at island locations
- Floor to ceiling models are also available
- One side of enclosure is removable for access to utilities
- Painted steel construction to match the cabinets
- Umbilicals are also available in Architectural Resin, Plastic Laminate, and Stainless Steel


### Island Bench Umbilicals

Note: A verified ceiling height is required to determine the correct umbilical height



Model	Width	Depth	Height	Model	Width	Depth	Height
UMI180665	18"	6"	65.25"	UMI240665	24"	6"	65.25"
UMI180669	18"	6"	69"	UMI240669	24"	6"	69"
UMI180671	18"	6"	71.25"	UMI240671	24"	6"	71.25"
UMI180675	18"	6"	75"	UMI240675	24"	6"	75"
UMI180677	18"	6"	77.25"	UMI240677	24"	6"	77.25"
UMI180681	18"	6"	81"	UMI240681	24"	6"	81"
UMI180683	18"	6"	83.25"	UMI240683	24"	6"	83.25"
UMI180689	18"	6"	89.25"	UMI240689	24"	6"	89.25"
UMI180693	18"	6"	93"	UMI240693	24"	6"	93"
UMI180695	18"	6"	95.25"	UMI240695	24"	6"	95.25"
UMI180865	18"	8"	65.25"	UMI240865	24"	8"	65.25"
UMI180869	18"	8"	69"	UMI240869	24"	8"	69"
UMI180871	18"	8"	71.25"	UMI240871	24"	8"	71.25"
UMI180875	18"	8"	75"	UMI240875	24"	8"	75"
UMI180877	18"	8"	77.25"	UMI240877	24"	8"	77.25"
UMI180881	18"	8"	81"	UMI240881	24"	8"	81"
UMI180883	18"	8"	83.25"	UMI240883	24"	8"	83.25"
UMI180889	18"	8"	89.25"	UMI240889	24"	8"	89.25"
UMI180893	18"	8"	93"	UMI240893	24"	8"	93"
UMI180895	18"	8"	95.25"	UMI240895	24"	8"	95.25"
UMI181065	18"	10"	65.25"	UMI241065	24"	10"	65.25"
UMI181069	18"	10"	69"	UMI241069	24"	10"	69"
UMI181071	18"	10"	71.25"	UMI241071	24"	10"	71.25"
UMI181075	18"	10"	75"	UMI241075	24"	10"	75"
UMI181077	18"	10"	77.25"	UMI241077	24"	10"	77.25"
UMI181081	18"	10"	81"	UMI241081	24"	10"	81"
UMI181083	18"	10"	83.25"	UMI241083	24"	10"	83.25"
UMI181089	18"	10"	89.25"	UMI241089	24"	10"	89.25"
UMI181093	18"	10"	93"	UMI241093	24"	10"	93"
UMI181095	18"	10"	95.25"	UMI241095	24"	10"	95.25"

### Wall Bench Umbilicals



Same as above except substitute "W" for the "I" in the part number prefix (UMI)

Example:

Model	Width	Depth	Height
UMW180665	18"	6"	65.25"



# WALL CABINETS

## Construction Details

### HLS Wall Cabinets





Hanson wall cabinets are available in a variety of styles and sizes designed to meet your specific applications. They are designed for maximum storage capacity, and access.


- 24", 30", 36" and 48" high wall cabinets
- Choice of 13" or 16" deep cabinets
- Available in solid panel, framed or frameless glass doors and open style wall cabinets
- Hinged or sliding door options



### HLS Features and Benefits


- Ⓐ Cabinet body design accepts multiple styles of doors.
- Ⓑ 18 ga. double door fronts, with sound deadening full channeled edges. Glass door applications feature 1/4" thick annealed clear glass.
- Ⓒ 18 ga. cabinet bodies with flush seamless welded face.
- Ⓓ 20 ga. adjustable shelves are included in open and cupboard cabinets.
- Ⓔ Stainless steel 5 knuckle hinges on cabinet doors that close against neoprene rubber bumpers.
- Ⓕ HLS cabinets have electrostatically applied, baked on powder coated epoxy paint finish for superior performance.
- Ⓖ Front shelf lips are included on all shelves.






## The HLS Advantages

- Inset steel - flush face construction
- Cabinets are 24" high - width and depth as shown below
- Furnished with 1 adjustable shelf




### Open Cabinets

Model	Depth	Width
G18EL	13"	18-3/4"
G18EL16	16"	18-3/4"
G24EL	13"	24-3/4"
G24EL16	16"	24-3/4"
G30EL	13"	30"
G30EL16	16"	30"
G35EL	13"	35"
G35EL16	16"	35"
G42EL	13"	42"
G42EL16	16"	42"
G47EL	13"	47"
G47EL16	16"	47"



### Hinged Solid Doors


Model	Depth	Width
A18EL	13"	18-3/4"
A18EL16	16"	18-3/4"
A24EL	13"	24-3/4"
A24EL16	16"	24-3/4"



### Left Hinged Solid Doors


Model	Depth	Width
A18ELL	13"	18-3/4"
A18ELL16	16"	18-3/4"
A24ELL	13"	24-3/4"
A24ELL16	16"	24-3/4"

L = Left Hinged Door




### Hinged Solid Doors

Model	Depth	Width
B30EL	13"	30"
B30EL16	16"	30"
B35EL	13"	35"
B35EL16	16"	35"
B42EL	13"	42"
B42EL16	16"	42"
B47EL	13"	47"
B47EL16	16"	47"




### Framed Glass Hinged Doors

Model	Depth	Width
D18EL	13"	18-3/4"
D18EL16	16"	18-3/4"
D24EL	13"	24-3/4"
D24EL16	16"	24-3/4"



### Sliding Solid Doors

Model	Depth	Width
C30EL	13"	30"
C30EL16	16"	30"
C35EL	13"	35"
C35EL16	16"	35"
C42EL	13"	42"
C42EL16	16"	42"
C47EL	13"	47"
C47EL16	16"	47"




### Framed Glass Left Hinged Doors

Model	Depth	Width
D18ELL	13"	18-3/4"
D18ELL16	16"	18-3/4"
D24ELL	13"	24-3/4"
D24ELL16	16"	24-3/4"

L = Left Hinged Door

# WALL CABINETS


Extra Low • 24"



## The HLS Advantages


- Inset steel - flush face construction
- Cabinets are 24" high - width and depth as shown below
- Furnished with 1 adjustable shelf

### Framed Glass Hinged Doors




Model	Depth	Width
E30EL	13"	30"
E30EL16	16"	30"
E35EL	13"	35"
E35EL16	16"	35"
E42EL	13"	42"
E42EL16	16"	42"
E47EL	13"	47"
E47EL16	16"	47"

### Frameless Glass Sliding Doors




Model	Depth	Width
FAL30EL	13"	30"
FAL30EL16	16"	30"
FAL35EL	13"	35"
FAL35EL16	16"	35"
FAL42EL	13"	42"
FAL42EL16	16"	42"
FAL47EL	13"	47"
FAL47EL16	16"	47"

### Framed Glass Sliding Doors




Model	Depth	Width
F30EL	13"	30"
F30EL16	16"	30"
F35EL	13"	35"
F35EL16	16"	35"
F42EL	13"	42"
F42EL16	16"	42"
F47EL	13"	47"
F47EL16	16"	47"



## The HLS Advantages

- Inset steel - flush face construction
- Cabinets are 30" high - width and depth as shown below
- Furnished with 2 adjustable shelf

### Open Cabinets




Model	Depth	Width
G18L	13"	18-3/4"
G18L16	16"	18-3/4"
G24L	13"	24-3/4"
G24L16	16"	24-3/4"
G30L	13"	30"
G30L16	16"	30"
G35L	13"	35"
G35L16	16"	35"
G42L	13"	42"
G42L16	16"	42"
G47L	13"	47"
G47L16	16"	47"

### Hinged Solid Doors




Model	Depth	Width
A18L	13"	18-3/4"
A18L16	16"	18-3/4"
A24L	13"	24-3/4"
A24L16	16"	24-3/4"

### Left Hinged Solid Doors




Model	Depth	Width
A18LL	13"	18-3/4"
A18LL16	16"	18-3/4"
A24LL	13"	24-3/4"
A24LL16	16"	24-3/4"

### Hinged Solid Doors




Model	Depth	Width
B30L	13"	30"
B30L16	16"	30"
B35L	13"	35"
B35L16	16"	35"
B42L	13"	42"
B42L16	16"	42"
B47L	13"	47"
B47L16	16"	47"

### Framed Glass Hinged Doors




Model	Depth	Width
D18L	13"	18-3/4"
D18L16	16"	18-3/4"
D24L	13"	24-3/4"
D24L16	16"	24-3/4"

### Sliding Solid Doors



Model	Depth	Width
C30L	13"	30"
C30L16	16"	30"
C35L	13"	35"
C35L16	16"	35"
C42L	13"	42"
C42L16	16"	42"
C47L	13"	47"
C47L16	16"	47"


### Framed Glass Left Hinged Doors



Model	Depth	Width
D18LL	13"	18-3/4"
D18LL16	16"	18-3/4"
D24LL	13"	24-3/4"
D24LL16	16"	24-3/4"

# WALL CABINETS


Low • 30"



## The HLS Advantages


- Inset steel - flush face construction
- Cabinets are 30" high - width and depth as shown below
- Furnished with 2 adjustable shelf

### Framed Glass Hinged Doors




Model	Depth	Width
E30L	13"	30"
E30L16	16"	30"
E35L	13"	35"
E35L16	16"	35"
E42L	13"	42"
E42L16	16"	42"
E47L	13"	47"
E47L16	16"	47"

### Frameless Glass Sliding Doors




Model	Depth	Width
FAL30L	13"	30"
FAL30L16	16"	30"
FAL35L	13"	35"
FAL35L16	16"	35"
FAL42L	13"	42"
FAL42L16	16"	42"
FAL47L	13"	47"
FAL47L16	16"	47"

### Framed Glass Sliding Doors




Model	Depth	Width
F30L	13"	30"
F30L16	16"	30"
F35L	13"	35"
F35L16	16"	35"
F42L	13"	42"
F42L16	16"	42"
F47L	13"	47"
F47L16	16"	47"



## The HLS Advantages

- Inset steel - flush face construction
- Cabinets are 36" high - width and depth as shown below
- Furnished with 2 adjustable shelf

### Open Cabinets




Model	Depth	Width
G18M	13"	18-3/4"
G18M16	16"	18-3/4"
G24M	13"	24-3/4"
G24M16	16"	24-3/4"
G30M	13"	30"
G30M16	16"	30"
G35M	13"	35"
G35M16	16"	35"
G42M	13"	42"
G42M16	16"	42"
G47M	13"	47"
G47M16	16"	47"

### Hinged Solid Doors




Model	Depth	Width
A18M	13"	18-3/4"
A18M16	16"	18-3/4"
A24M	13"	24-3/4"
A24M16	16"	24-3/4"

### Left Hinged Solid Doors




Model	Depth	Width
A18ML	13"	18-3/4"
A18ML16	16"	18-3/4"
A24ML	13"	24-3/4"
A24ML16	16"	24-3/4"

### Hinged Solid Doors




Model	Depth	Width
B30M	13"	30"
B30M16	16"	30"
B35M	13"	35"
B35M16	16"	35"
B42M	13"	42"
B42M16	16"	42"
B47M	13"	47"
B47M16	16"	47"




### Framed Glass Hinged Doors

Model	Depth	Width
D18M	13"	18-3/4"
D18M16	16"	18-3/4"
D24M	13"	24-3/4"
D24M16	16"	24-3/4"

### Sliding Solid Doors



Model	Depth	Width
C30M	13"	30"
C30M16	16"	30"
C35M	13"	35"
C35M16	16"	35"
C42M	13"	42"
C42M16	16"	42"
C47M	13"	47"
C47M16	16"	47"




### Framed Glass Left Hinged Doors

Model	Depth	Width
D18ML	13"	18-3/4"
D18ML16	16"	18-3/4"
D24ML	13"	24-3/4"
D24ML16	16"	24-3/4"

# WALL CABINETS


Medium • 36"




## The HLS Advantages

- Inset steel - flush face construction
- Cabinets are 36" high - width and depth as shown below
- Furnished with 2 adjustable shelf

### Framed Glass Hinged Doors




Model	Depth	Width
E30M	13"	30"
E30M16	16"	30"
E35M	13"	35"
E35M16	16"	35"
E42M	13"	42"
E42M16	16"	42"
E47M	13"	47"
E47M16	16"	47"



### Frameless Glass Sliding Doors

Model	Depth	Width
FAL30M	13"	30"
FAL30M16	16"	30"
FAL35M	13"	35"
FAL35M16	16"	35"
FAL42M	13"	42"
FAL42M16	16"	42"
FAL47M	13"	47"
FAL47M16	16"	47"


### Framed Glass Sliding Doors



Model	Depth	Width
F30M	13"	30"
F30M16	16"	30"
F35M	13"	35"
F35M16	16"	35"
F42M	13"	42"
F42M16	16"	42"
F47M	13"	47"
F47M16	16"	47"

High • 48"


[hansonlab.com](http://hansonlab.com) • 805-498-3121



## The HLS Advantages


- Inset steel - flush face construction
- Cabinets are 48" high - width and depth as shown below
- Furnished with 3 adjustable shelf

### Open Cabinets




Model	Depth	Width
G18H	13"	18-3/4"
G18H16	16"	18-3/4"
G24H	13"	24-3/4"
G24H16	16"	24-3/4"
G30H	13"	30"
G30H16	16"	30"
G35H	13"	35"
G35H16	16"	35"
G42H	13"	42"
G42H16	16"	42"
G47H	13"	47"
G47H16	16"	47"

### Hinged Solid Doors




Model	Depth	Width
A18H	13"	18-3/4"
A18H16	16"	18-3/4"
A24H	13"	24-3/4"
A24H16	16"	24-3/4"

### Left Hinged Solid Doors




Model	Depth	Width
A18HL	13"	18-3/4"
A18HL16	16"	18-3/4"
A24HL	13"	24-3/4"
A24HL16	16"	24-3/4"

### Hinged Solid Doors



Model	Depth	Width
B30H	13"	30"
B30H16	16"	30"
B35H	13"	35"
B35H16	16"	35"
B42H	13"	42"
B42H16	16"	42"
B47H	13"	47"
B47H16	16"	47"


### Sliding Solid Doors



Model	Depth	Width
C30H	13"	30"
C30H16	16"	30"
C35H	13"	35"
C35H16	16"	35"
C42H	13"	42"
C42H16	16"	42"
C47H	13"	47"
C47H16	16"	47"

# WALL CABINETS


High • 48"



## The HLS Advantages


- Inset steel - flush face construction
- Cabinets are 48" high - width and depth as shown below
- Furnished with 3 adjustable shelf

### Framed Glass Hinged Doors




Model	Depth	Width
D18H	13"	18-3/4"
D18H16	16"	18-3/4"
D24H	13"	24-3/4"
D24H16	16"	24-3/4"

### Framed Glass Left Hinged Doors




Model	Depth	Width
D18HL	13"	18-3/4"
D18HL16	16"	18-3/4"
D24HL	13"	24-3/4"
D24HL16	16"	24-3/4"

### Framed Glass Hinged Doors




Model	Depth	Width
E30H	13"	30"
E30H16	16"	30"
E35H	13"	35"
E35H16	16"	35"
E42H	13"	42"
E42H16	16"	42"
E47H	13"	47"
E47H16	16"	47"

### Framed Glass Sliding Doors



Model	Depth	Width
F30H	13"	30"
F30H16	16"	30"
F35H	13"	35"
F35H16	16"	35"
F42H	13"	42"
F42H16	16"	42"
F47H	13"	47"
F47H16	16"	47"


### Frameless Sliding Doors



Model	Depth	Width
FAL30H	13"	30"
FAL30H16	16"	30"
FAL35H	13"	35"
FAL35H16	16"	35"
FAL42H	13"	42"
FAL42H16	16"	42"
FAL47H	13"	47"
FAL47H16	16"	47"

## Wall Cabinet Front Filler Panel


- For use between 2 cabinets or between cabinet & wall
- Specify T if the panel needs to be field trimmed with an SA1WC double-slip angle



Model	Width	Height	Depth
WCF1EL13	1"	24"	13"
WCF1L13	1"	30"	13"
WCF1M13	1"	36"	13"
WCF1H13	1"	48"	13"
WCF3EL13	3"	24"	13"
WCF3L13	3"	30"	13"
WCF3M13	3"	36"	13"
WCF3H13	3"	48"	13"
WCF6EL13	6"	24"	13"
WCF6L13	6"	30"	13"
WCF6M13	6"	36"	13"
WCF6H13	6"	48"	13"
WCF12EL13	12"	24"	13"
WCF12L13	12"	30"	13"
WCF12M13	12"	36"	13"
WCF12H13	12"	48"	13"
WCF1EL16	1"	24"	16"
WCF1L16	1"	30"	16"
WCF1M16	1"	36"	16"
WCF1H16	1"	48"	16"
WCF3EL16	3"	24"	16"
WCF3L16	3"	30"	16"
WCF3M16	3"	36"	16"
WCF3H16	3"	48"	16"
WCF6EL16	6"	24"	16"
WCF6L16	6"	30"	16"
WCF6M16	6"	36"	16"
WCF6H16	6"	48"	16"
WCF12EL16	12"	24"	16"
WCF12L16	12"	30"	16"
WCF12M16	12"	36"	16"
WCF12H16	12"	48"	16"

## Wall Cabinet Corner Filler Panel


- Used to close top and bottom where 2 upper cabinets are joined at right angles
- Supplied in sets of two



Model	Width	Depth
WCCF13	13"	13"
WCCF16	16"	16"
WCCF1316	13"	16"

## Wall Cabinet Slip Angle

- Field trimmable for all height wall cabinets



Model	Width	Height
SA1WC24	1"	24"
SA1WC30	1"	30"
SA1WC36	1"	36"
SA1WC48	1"	48"


Includes a single angle.

# WALL CABINET ACCESSORIES

## Ceiling Enclosures & Sloping Tops

### Wall Cabinet Ceiling Enclosures


- Used to enclose the space between the top of the wall cabinet and the ceiling
- Angles for installation and trimmable scribes are provided
- A verified ceiling height is required to determine the correct enclosure height



Model	Width	Depth	Height
CEW1813-18	18.75"	13"	18"
CEW2413-18	24.75"	13"	18"
CEW3013-18	30"	13"	18"
CEW3513-18	35"	13"	18"
CEW4213-18	42"	13"	18"
CEW4713-18	47"	13"	18"
CEW1816-18	18.75"	16"	18"
CEW2416-18	24.75"	16"	18"
CEW3016-18	30"	16"	18"
CEW3516-18	35"	16"	18"
CEW4216-18	42"	16"	18"
CEW4716-18	47"	16"	18"
CEW1813-24	18.75"	13"	24"
CEW2413-24	24.75"	13"	24"
CEW3013-24	30"	13"	24"
CEW3513-24	35"	13"	24"
CEW4213-24	42"	13"	24"
CEW4713-24	47"	13"	24"
CEW1816-24	18.75"	16"	24"
CEW2416-24	24.75"	16"	24"
CEW3016-24	30"	16"	24"
CEW3516-24	35"	16"	24"
CEW4216-24	42"	16"	24"
CEW4716-24	47"	16"	24"
CEW1813-30	18.75"	13"	30"
CEW2413-30	24.75"	13"	30"
CEW3013-30	30"	13"	30"
CEW3513-30	35"	13"	30"
CEW4213-30	42"	13"	30"
CEW4713-30	47"	13"	30"
CEW1816-30	18.75"	16"	30"
CEW2416-30	24.75"	16"	30"
CEW3016-30	30"	16"	30"
CEW3516-30	35"	16"	30"
CEW4216-30	42"	16"	30"
CEW4716-30	47"	16"	30"

### Sloping Tops


- Used to enclose the space between the top of the wall cabinet and the ceiling



Model	Width	Depth	Height
A4810-18	18.75"	13"	8.344"
A4810-24	24.75"	13"	8.344"
A4810-30	30"	13"	8.344"
A4810-35	35"	13"	8.344"
A4810-42	42"	13"	8.344"
A4810-47	47"	13"	8.344"
A4820-18	18.75"	16"	9.910"
A4820-24	24.75"	16"	9.910"
A4820-30	30"	16"	9.910"
A4820-35	35"	16"	9.910"
A4820-42	42"	16"	9.910"
A4820-47	47"	16"	9.910"

### Corner Sloping Top

- Used in corners with sloping tops and where wall cabinets are at right angles





Model	Depth	Height
CSL13	13"	8.344"
CSL16	16"	9.910"

### HLS Tall Cabinets


Hanson tall cabinets are available in a variety of styles and sizes designed to meet your specific applications. They are designed for maximum storage capacity, and access.

- Tall cabinets are 84" high
- Choice of 18" or 23" deep cabinets
- Available in solid panel or framed glass in both hinged and sliding doors
- Instrument and wardrobe cabinets are also available




### HLS Features and Benefits

- Ⓐ Cabinet body design accepts multiple styles of doors.
- Ⓑ 18 ga. double door fronts, with sound deadening full channeled edges. Glass door applications feature 1/4" thick annealed clear glass.
- Ⓒ 18 ga. cabinet bodies with flush seamless welded face.
- Ⓓ 5 each 20 ga. adjustable shelves are included in open and door cabinets. Shelves have welded hat-channels in cabinets 35" and over.
- Ⓔ Stainless steel 5 knuckle hinges on cabinet doors that close against neoprene rubber bumpers.
- Ⓕ HLS cabinets have electrostatically applied, baked on powder coated epoxy paint finish for superior performance.
- Ⓖ 14 ga. gusset plates with hex-head leveling bolts at each corner.
- Ⓗ Front shelf lips are included on all shelves.




# TALL STORAGE CABINETS

Height • 84"




## Solid Door - Tall Cabinets

- Inset steel - flush face construction
- Tall storage cabinets are 84" high - width and depth as shown below
- Tall cabinets have 5 adjustable shelves




### Hinged Doors

Model	Depth	Width
H18HA	18"	18-3/4"
H18HB	23"	18-3/4"
H24HA	18"	24-3/4"
H24HB	23"	24-3/4"




### Left Hinged Doors

Model	Depth	Width
H18HAL	18"	18-3/4"
H18HBL	23"	18-3/4"
H24HAL	18"	24-3/4"
H24HBL	23"	24-3/4"



### Hinged Doors

Model	Depth	Width
I30HA	18"	30"
I30HB	23"	30"
I35HA	18"	35"
I35HB	23"	35"
I42HA	18"	42"
I42HB	23"	42"
I47HA	18"	47"
I47HB	23"	47"




### Sliding Doors

Model	Depth	Width
J30HA	18"	30"
J30HB	23"	30"
J35HA	18"	35"
J35HB	23"	35"
J42HA	18"	42"
J42HB	23"	42"
J47HA	18"	47"
J47HB	23"	47"


Height • 84"

[hansonlab.com](http://hansonlab.com) • 805-498-3121




## Framed Glass Door - Tall Cabinets

- Inset steel - flush face construction
- Tall storage cabinets are 84" high - width and depth as shown below
- Tall cabinets have 5 adjustable shelves




### Hinged Doors

Model	Depth	Width
K18HA	18"	18-3/4"
K18HB	23"	18-3/4"
K24HA	18"	24-3/4"
K24HB	23"	24-3/4"




### Left Hinged Doors

Model	Depth	Width
K18HAL	18"	18-3/4"
K18HBL	23"	18-3/4"
K24HAL	18"	24-3/4"
K24HBL	23"	24-3/4"



### Hinged Doors

Model	Depth	Width
L30HA	18"	30"
L30HB	23"	30"
L35HA	18"	35"
L35HB	23"	35"
L42HA	18"	42"
L42HB	23"	42"
L47HA	18"	47"
L47HB	23"	47"




### Sliding Doors

Model	Depth	Width
M30HA	18"	30"
M30HB	23"	30"
M35HA	18"	35"
M35HB	23"	35"
M42HA	18"	42"
M42HB	23"	42"
M47HA	18"	47"
M47HB	23"	47"


# TALL STORAGE CABINETS

Height • 84"



## Open - Tall Cabinets


- Inset steel - flush face construction.
- Tall storage cabinets are 84" high - width and depth as shown below
- Tall cabinets have 5 adjustable shelves



Model	Depth	Width
N18HA	18"	18-3/4"
N18HB	23"	18-3/4"
N24HA	18"	24-3/4"
N24HB	23"	24-3/4"
N30HA	18"	30"
N30HB	23"	30"
N35HA	18"	35"
N35HB	23"	35"
N42HA	18"	42"
N42HB	23"	42"
N47HA	18"	47"
N47HB	23"	47"


Height • 84"

[hansonlab.com](http://hansonlab.com) • 805-498-3121




## Combination - Tall Cabinet

- Inset steel - flush face construction
- Tall storage cabinets are 84" high - width and depth as shown below
- Solid hinged door at bottom
- Top cabinet has (3) adjustable shelves, bottom cabinet has (1) adjustable shelf



### Framed Glass Hinged Door at Top

Model	Depth	Width
O18HA	18"	18-3/4"
O18HB	23"	18-3/4"
O24HA	18"	24-3/4"
O24HB	23"	24-3/4"




### Left Framed Glass Hinged Door at Top

Model	Depth	Width
O18HAL	18"	18-3/4"
O18HBL	23"	18-3/4"
O24HAL	18"	24-3/4"
O24HBL	23"	24-3/4"



### Framed Glass Hinged Doors at Top

Model	Depth	Width
O30HA	18"	30"
O30HB	23"	30"
O35HA	18"	35"
O35HB	23"	35"
O42HA	18"	42"
O42HB	23"	42"
O47HA	18"	47"
O47HB	23"	47"




### Framed Glass Sliding Doors at Top

Model	Depth	Width
P30HA	18"	30"
P30HB	23"	30"
P35HA	18"	35"
P35HB	23"	35"
P42HA	18"	42"
P42HB	23"	42"
P47HA	18"	47"
P47HB	23"	47"


# TALL STORAGE CABINETS

Height • 84"




## Wardrobe - Tall Cabinets

- Inset steel - flush face construction
- Tall storage cabinets are 84" high - width and depth as shown below
- Top Shelf
- Clothes hanger rod
- Hinged solid door




### Hinged Doors

Model	Depth	Width
Q18HB	23"	18-3/4"
Q24HB	23"	24-3/4"



### Left Hinged Doors


Model	Depth	Width
Q18HBL	23"	18-3/4"
Q24HBL	23"	24-3/4"



### Hinged Doors

Model	Depth	Width
Q30HB	23"	30"
Q35HB	23"	35"
Q42HB	23"	42"
Q47HB	23"	47"

## Tall Cabinet Front / End Filler



- Used between two cabinets or between a cabinet and a wall
- Specify -T if the panel needs to be field trimmed with an SA1TC double-slip angle

### Front Filler

Model	Width	Height
TCF1	1"	84"
TCF3	3"	84"
TCF6	6"	84"
TCF12	12"	84"

### End Filler


Model	Width	Height
TCEF1	1"	84"

### Note:

TCF1 & TCEF1 are not trimmable.

## Tall Cabinet Finished Backs


- Finished back of tall cabinets
- Finished back panels are attached at the factory
- Add -FB to the end of the tall cabinet model number
- Tall cabinets must be properly anchored to walls or floors to prevent tipping



Model  
-FB

## Double-Slip Angle


- Used with tall cabinet when front filler panels need to be field trimmed



Model	Width	Height
SA1TC	1"	84"

# TALL CABINET ACCESSORIES

## Ceiling Enclosures & Sloping Tops





### Tall Cabinet Ceiling Enclosures

- Used to enclose the space between the top of the wall cabinet and the ceiling
- Angles for installation and trimmable scribes are provided
- A verified ceiling height is required to determine the correct enclosure height

### Sloping Tops


- Add -SL to the end of the tall cabinet model number



Model	Width	Depth	Height
A6810-18	18.75"	18"	11.064"
A6810-24	24.75"	18"	11.064"
A6810-30	30"	18"	11.064"
A6810-35	35"	18"	11.064"
A6810-42	42"	18"	11.064"
A6810-47	47"	18"	11.064"
A6820-18	18.75"	23"	13.951"
A6820-24	24.75"	23"	13.951"
A6820-30	30"	23"	13.951"
A6820-35	35"	23"	13.951"
A6820-42	42"	23"	13.951"
A6820-47	47"	23"	13.951"

### Corner Sloping Top

- Used in corners with sloping tops and where tall cabinets are at right angles



Model	Depth	Height
CSL18	18"	11.064"
CSL23	23"	13.951"

Model	Width	Depth	Height
CET1818-24	18.75"	18"	24"
CET1823-24	18.75"	23"	24"
CET2418-24	24.75"	18"	24"
CET2423-24	24.75"	23"	24"
CET3018-24	30"	18"	24"
CET3023-24	30"	23"	24"
CET3518-24	35"	18"	24"
CET3523-24	35"	23"	24"
CET4218-24	42"	18"	24"
CET4223-24	42"	23"	24"
CET4718-24	47"	18"	24"
CET4723-24	47"	23"	24"
CET1818-30	18.75"	18"	30"
CET1823-30	18.75"	23"	30"
CET2418-30	24.75"	18"	30"
CET2423-30	24.75"	23"	30"
CET3018-30	30"	18"	30"
CET3023-30	30"	23"	30"
CET3518-30	35"	18"	30"
CET3523-30	35"	23"	30"
CET4218-30	42"	18"	30"
CET4223-30	42"	23"	30"
CET4718-30	47"	18"	30"
CET4723-30	47"	23"	30"
CET1818-36	18.75"	18"	36"
CET1823-36	18.75"	23"	36"
CET2418-36	24.75"	18"	36"
CET2423-36	24.75"	23"	36"
CET3018-36	30"	18"	36"
CET3023-36	30"	23"	36"
CET3518-36	35"	18"	36"
CET3523-36	35"	23"	36"
CET4218-36	42"	18"	36"
CET4223-36	42"	23"	36"
CET4718-36	47"	18"	36"
CET4723-36	47"	23"	36"
CET1818-48	18.75"	18"	48"
CET1823-48	18.75"	23"	48"
CET2418-48	24.75"	18"	48"
CET2423-48	24.75"	23"	48"
CET3018-48	30"	18"	48"
CET3023-48	30"	23"	48"
CET3518-48	35"	18"	48"
CET3523-48	35"	23"	48"
CET4218-48	42"	18"	48"
CET4223-48	42"	23"	48"
CET4718-48	47"	18"	48"
CET4723-48	47"	23"	48"