


M3 Series Workstation (Single-Sided) Specification Worksheet

M4 Series Workstation (Double-Sided) Specification Worksheet

Application notes:

- Max 5 duplexes per power cord.
- Max 2 power cords per post.
- Hose and Cables typically extend out of the top of the post
- Power and Gas **cannot** be shared in same post
- Gas and Data **can** share post.
- Power and Data **can** share post.



Nominal Bench Width

- ☐ 36"
- ☐ 48"
- ☐ 60"
- ☐ 72"

Nominal Bench Depth



- ☐ 30"
- ☐ 36"
- ☐ 60" (Double-Sided)

Optional



- ☐ Suspended cabinet rail kit
- ☐ Work surface, 1" thick Black Phenolic
- ☐ Work surface other:

Plumbing Services Pre-Plumbed with quick connects (choose one plate style per post)



M3 Series Workstation (Single-Sided)



M4 Series Workstation (Double-Sided)





Below-Counter Data, (default is None; max 1 total per post)



Electrical Services (choose one plate style per post)
See page 4 for cord/plug information


Left Post


None	
None	
None	
None	


None		None
None		None
None		None
None		None

None	
-----------------------------------	---


Right Post


	None
	None
	None
	None

None		None
None		None
None		None
None		None

	None
--	-----------------------------------

Below-Counter Power, (default is None; max 1 total per post)


None	
-----------------------------------	---

	None
---	-----------------------------------


Plumbing Services

Pre-plumbed with quick connects


Needle Valve - Black Knob (default valve)



Needle Valve - "ColorTech" Knob



Ball Valve - Lever Handle



Electrical Services -

Typical Electrical Plug Detail

Each cable is designed to support a single circuit and comes standard with a NEMA L5-20P locking plug (one cord per circuit). Assemblies include an additional 10 feet of cord slack for flexible installation.

For applications requiring two circuits with a shared neutral, a single cable is equipped with a NEMA L14-20P locking plug. These assemblies also include 10 feet of additional cord slack.

Single-Circuit 20A, L5-20P, 1 cord

NEMA L5-20P
20A-125V (default plug)



NEMA L14-20P
20A-125/250V



NEMA L6-30P
30A-250V



Decorator Style Duplex Color Options



Brown
Gray
Orange
Red
White

120VAC, 20Amp Duplex



White

120VAC, 20Amp GFCI Duplex



White

5V output USB, 4-port