

Which M-Series® bench is right for your lab?

Configurable lab workstations designed for a range of layouts, equipment loads, and workflow requirements.

						
Model	M1	M2	M3	M4	M5	M6
Configuration	Single-sided	Double-sided	Single-sided	Double-sided	Single-sided	Single-sided
Construction	14-ga. 2" x 3" system	14-ga. 2" x 3" system	11-ga. 2" x 5" system	11-ga. 2" x 5" system	11-ga. 2" x 5" system	11-ga. 2" x 5" system
Duty Level	Mid-level	Mid-level	Heavy-duty	Heavy-duty	Heavy-duty	Heavy-duty
Mobility	Moveable with leg glides*	Moveable with leg glides*	Moveable with leg glides	Moveable with leg glides	Mobile	Mobile
Widths	36, 48, 60, 72 in.	36, 48, 60, 72 in.	36, 48, 60, 72 in.	36, 48, 60, 72 in.	36, 48, 60, 72 in.	36, 48, 60, 72 in.
Depths	30, 36 in.	60 in.	30, 36 in.	60 in.	30, 36, 42 in.	30, 36, 42 in.
Work Surface Height	Adjustable, 30–36 in.	Adjustable, 30–36 in.	Adjustable, 30–36 in.	Adjustable, 30–36 in.	Fixed, 30, 32, or 36 in.	Adjustable, 30–36 in.
Load Capacity	1,200 lb.	2,400 lb. total (1,200 lb. per side)	1,200 lb.	2,400 lb. total (1,200 lb. per side)	1,200 lb.	1,200 lb.
Power	Up to 5 duplexes per post; 1 cord per post	Up to 5 duplexes per post; 1 cord per post	Up to 9 duplexes per post; 4 cords per post	Up to 9 duplexes per post; 4 cords per post	Up to 9 duplexes per post; 4 cords per post	Up to 9 duplexes per post; 4 cords per post

* Optional caster configuration available for light-duty applications.

Build Your M-Series® Workstation

Integrated Power, Data & Utility Options

Integrated Power

Keep equipment connections accessible and organized.



Specialty Gas

Integrate pre-plumbed gas services into select workstation configurations.

Data & Connectivity

Configure data options based on model and project needs.

Durable Powder-Coat Finish Options



Petal White

Light Gray

Dark Gray

Shadow


Mountain Blue

Blue

Black

Custom color options may be available upon request. Additional lead time and charges may apply.

Configure Your Workstation



Work Surfaces

Phenolic and epoxy options available based on application requirements.



Shelving

Multiple configurations support equipment, supplies, and workflow needs.



Storage

Mobile and suspended cabinet options vary by model.



Accessories

Add power strips, bin box mounts and project-specific accessories.

Need help selecting the right M-Series® configuration?

Our lab planning specialists can help evaluate your layout, equipment loads, utility requirements, and workflow needs.

Connect with a Lab Planning Specialist →

www.hansonlab.com